
# TRF : Time-Resolved Fluoremetry

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
	Eu-chelate	Fluorescein
Excitation (nm)	340	490
Emission max. (nm)	613	515
Stokes' shift (nm)	273	25
Decay time (ns)	730,000	3
Emission width	Narrow	Broad
Background interference	Low	High
Self quenching	No	Strong
# of repeated measurements	High	Low
Applicability of multi-label	Good	Limited

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


DELFLIA：dissociation enhanced lanthanide FIA

LANCE：time-resolved fluorescent resonance energy transfer assay

TruPoint：time-resolved fluorescence quench assay

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##### 1. Immunoassays

- Madersbacher S, Shu-Chen T, Schwarz S et al. "Time-resolved immunofluorometry and other frequently used immunoassay types for follicle-stimulating hormone compared by using identical monoclonal antibodies" *Clin Chem* 1993; 39: 1435-9.

2. Kinase assays
  - Gaarde WA, etc., “Development of a Nonradioactive, Time-Resolved Fluorescence Assay for the Measurement of Jun N-terminal Kinase Activity” *J Biomol Screen* 1997 Aug **2:4** 213-23
  
3. Receptor-ligand binding assays
  - Inglese J, etc., “Chemokine receptor-ligand interactions measured using time-resolved fluorescence” *Biochemistry* 1998 Feb 24 **37:8** 2372-7
  - Moore KJ, etc., “A Homogenous 384-Well High Throughput Screening for Novel Tumor Necrosis Factor Receptor: Ligand Interactions Using Time Resolved Energy Transfer” *J Biomol Screen* 1999 Aug **4:4** 205-14
  - Makishima M, etc., “Identification of a nuclear receptor for bile acids” *Science* 1999 May 21 284:5418 1362-5
  
4. GPCR function assays: ex. cAMP, GTP binding
  - Valenzano, K.J., Miller, W., Kravitz J. N., Samama, P., Fitzpatrick, D., and Seeley, K. “Development of a fluorescent ligand-binding assay using the AcroWell filter plate” *J Biomol Screen* 2000, **5**:455-461
  
5. DNA hybridization
  - Sjooroo M, etc., “Triple-label hybridization assay for type-1 diabetes-related HLA alleles” *Biotechniques* 1995 May 18:5 870-7
  
6. Cytokine assay
  - Stenroos K, etc., “Homogeneous time-resolved IL-2-IL-2R alpha assay using fluorescence resonance energy transfer” *Cytokine* 1998 Jul 10:7 495-9
  
7. Biodistribution
  - Neville ME, etc., “A comparison of biodistribution of liposomal and soluble IL-2 by a new method based on time-resolved fluorometry of europium” *Cytokine* 2000 Nov 12:11 1702-11
  
8. Molecular interaction
  - Hamy F, etc., “Merged Screening for Human Immunodeficiency Virus Tat and Rev Inhibitors” *J Biomol Screen* 2001 Jun **6:3** 179-87
  
9. Enzyme assays
  - Boisclair MD, etc., “Development of a ubiquitin transfer assay for high throughput screening by fluorescence resonance energy transfer” *J Biomol Screen* 2000 Oct 5:5 319-28