


HUMAN HEALTH

ENVIRONMENTAL HEALTH

THE GPCR

PORTFOLIO THAT LEAVES
NOTHING OUT



GPCR
APPLICATION GUIDE

TABLE OF CONTENTS

Select Your Total GPCR Solution	4
GPCR Cell Lines	
Stable Recombinant Cell Lines	6
Calcium Detection Photoprotein Cell Lines	10
Ready-to-Use Frozen Cells	12
Cellular Detection Assays and Reagents	
cAMP Cellular Screening	14
Phospho-ERK Assays	18
Reporter Gene Assays	20
Label-free Validation	22
Binding Assays and Reagents	
GPCR Membrane Preparations	24
Ligand Receptor Binding	24
GTP Binding Assays	26
High Content Screening	28
Detection Instruments and Automation	
Instrumentation	32
Multilabel Plate Readers	33
Ultra HTS Microplate Imagers	34
Liquid Scintillation and Luminescence Counters	35
Automatic Gamma Counters	35
Automated Workstations	36
High Content Imagers	37
Microplates	38
Microplate Quick Reference Guide by Detection Method	40
Cell Line Cross Reference Guide	42



BIGGER GPCR SOLUTIONS

Get more from your research


With much drug discovery focused on GPCR targets, the need for reliable, proven and versatile assays, reagents and instrumentation to understand these complex signaling mechanisms has never been greater.

The PerkinElmer team understands the significant challenges of identifying biologically and pharmacologically relevant pathways. Scientist to scientist, we have the depth of experience in both cellular and biochemical platforms to support your work every step of the way.

We partner with you to quickly optimize your research, customizing solutions for specific target research.

From a portfolio that leaves nothing out

PerkinElmer offers the widest selection of proven GPCR target solutions – beyond the ligand GPCR event to the actual assessment of the cellular-activated signal transduction cascade. Enabling you to identify, validate, screen and profile your leads is our priority.



SELECT YOUR TOTAL GPCR SOLUTION

PerkinElmer is unequalled in providing an integrated solution that accelerates your research and advances drug discovery.

Right from the start, we help you determine the best strategy to identify and monitor:

- Effects of GPCR activation: changes in intracellular calcium, cAMP, phospho-ERK, inositol phosphate
- Reporter gene expression
- Receptor ligand binding
- GTP binding

Measurement	Intracellular Ca ²⁺		Inositol Phosphates	Intracellular cAMP		ERK Phosphorylation
Pathway	Phospholipase C		Phospholipase C	Adenylyl cyclase (AC)		MAP kinase
Assay Type	Functional	Image-based single-cell readout	Functional	Functional		Functional
G-protein Coupling	Optimal for G _s -coupled receptors, as well as for many G _i -coupled receptors. Can be used with G _i - and G _s -coupled receptors in conjunction with promiscuous G protein when needed.		Optimal for G _s -coupled receptors, as well as for many G _i -coupled receptors. Can be used with G _i - and G _s -coupled receptors in conjunction with promiscuous G protein when needed.	G _i (inhibitory), G _s (stimulatory)		G _s , G _i , G _q (in some instances)
Technologies	Fluorescent dyes (many AequoScreen and ValiScreen cell lines are validated for this assay). Photoproteins for flash luminescence assays: AequoScreen, PhotoScreen®, AequoZen, microplates.	Calcium flux using calcium fluorescent sensors, microplates	IP ₃ measurements: SPA beads, radiolabeled IP ₃ precursor, AequoScreen® cell lines, AequoZen, microplates	LANCE Ultra cAMP, AlphaScreen cAMP, cAMPZen frozen cells and ValiScreen cell lines validated for cAMP measurement, microplates	NEN SPA homogeneous radioimmunoassay, cAMP and adenylyl cyclase FlashPlate, microplates	AlphaScreen® SureFire® Phospho-ERK 1/2, AlphaScreen® SureFire® Total ERK, microplates
Advantages	Highly sensitive automated screening system to detect agonists and antagonists. Real-time kinetic readout.	Luminescence is extremely sensitive and cost-effective, enabling detection of partial agonists as well as allosteric modulators. Fluorescence is widely adopted and can be used when flash luminescence readers are not available.	Direct IP measurement, useful for confirmation of calcium flux hits using another type of assay	Highly sensitive automated screening system to detect agonists and antagonists. Detects elevations in basal cAMP due to constitutive activity and thus is able to detect inverse agonist activity. Non-radioactive.	Most sensitive ELISA format. Binding signal generation is a direct signal, not indirect like other assays. Less complex and better suited when therapeutic antibodies are developed.	Generic GPCR technology for G _s , G _i and G _q -coupled GPCRs. Excellent for receptors optimally coupled to calcium or cAMP pathway. Positive readout for G _i . Can be used with endogenously expressed receptors and stem cells.
Reader Requirements	Luminescence or fluorescence reader with injectors	High content screening readers with injectors	Radiometric detector	TR-FRET/Alpha-enabled reader	Radiometric detector	Alpha-enabled reader
PerkinElmer Readers	EnVision® Multilabel Plate Reader with Injectors, VICTOR™ X Multilabel Plate Reader with Dispenser, MicroBeta™ LumIJET™ Microplate Counter	Opera™ HCS System with Dispenser, Opera LX HCS System with Dispenser and UV Upgrade	ViewLux™ uHTS Imager, MicroBeta™ Scintillation & Luminescence Counter, TopCount™ Scintillation & Luminescence Counter, Tri-Carb® Liquid Scintillation Analyzer	ViewLux uHTS Imager, EnVision Multilabel Plate Reader with TRF Laser or Alpha Reader, VICTOR X4 or X5 Multilabel Plate Reader, EnSpire™ Multilabel Plate Reader	ViewLux uHTS Imager, MicroBeta™ Scintillation & Luminescence Counter, TopCount Scintillation & Luminescence Counter	EnVision Multilabel Plate Reader with Alpha Option, EnSpire Multilabel Plate Reader
Page Number	33, 34, 35	37	34, 35, 36	33, 34	34, 35	33

*Comparison of cAMP Assay Technologies for High Throughput Screening. Patricia Kasila and Harry Harney, PerkinElmer, Inc. Scientific Poster 4083.

Our state-of-the-art solutions portfolio is extensive and includes:

- Scalable cellular imaging, detection instrumentation and liquid handling systems for radiometric, fluorescence, luminescence or photometric detection
- Optimized assay platforms
- Versatile reagents – GPCR cell lines, frozen ready-to-use cells, radioligands, ligands and cloned receptor membranes

Innovation, combined with proven expertise and around-the-clock service, results in best-in-class solutions whatever the need. You can be confident that PerkinElmer specialists will help you find the most effective, and efficient, discovery strategy.

Reporter Gene Activity NFAT (Ca ²⁺ signaling), CRE (cAMP signaling) or Others	GTPγS Binding	GPCR Internalization		β-arrestin Recruitment	Receptor Ligand Binding	
Phospholipase C/AC/others	G-protein activation	–		–	–	
Functional	Functional	Functional, image-based readout, receptor endocytosis		Functional, image-based readout, clathrin-coated pit/vesicle formation	Biochemical	
Generic	G _i and some G _s , G _q	G _s , G _q , G ₁₂ , G ₁₃ as well as non-G protein-coupled GPCRs		G _s , G _q , G ₁₂ , G ₁₃ as well as non-G protein-coupled GPCRs Ideal for orphan GPCRs	Generic	
britelite plus, neolite, steadylite plus, nuclear translocation of NFAT using fluorescently labeled antibodies, microplates	SPA WGA-coated beads, NEN radionuclides, membrane preparations, radioligands, microplates	Homemade GFP-tagged receptor, GFP-tagged antibodies, microplates		Clathrin-coated pit/vesicle formation using GFP-tagged β-arrestin, microplates	NEN radioligands, SPA WGA-coated beads, WGA FlashPlate, membrane preparations, microplates	DELFI ligands
Highly sensitive, cost-effective GPCR screening system	Functional assay for GPCR activation, close to receptor activation	No artifacts due to varying cell number (transient transfections, pipetting errors)		Highly sensitive follow-up screening system to detect agonists and antagonists. No artifacts due to varying cell number (transient transfections, pipetting errors).	Saturation binding and competitive binding assays possible. Widely used technology for GPCR profiling and screening.	
Long stimulation time may reveal interesting pharmacologies	Allows detection of inverse agonism					
Luminescence reader	Radiometric detector	High content screening reader	Radiometric detector	High content screening reader	Radiometric detector	TRF-enabled reader
EnVision Multilabel Plate Reader, VICTOR X Multilabel Plate Reader, MicroBeta ² Scintillation & Luminescence Counter, TopCount Scintillation & Luminescence Counter, ViewLux uHTS Imager, TopCount Scintillation & Luminescence Counter, EnSpire Multilabel Plate Reader, Opera, Opera LX and Operetta [®] HCS systems	ViewLux uHTS Imager, MicroBeta ² Scintillation & Luminescence Counter, TopCount Scintillation & Luminescence Counter, Tri-Carb Liquid Scintillation Analyzer	Opera HCS System, Opera LX HCS, Operetta HCS Reader, UltraVIEW [®] VoX	TopCount Scintillation & Luminescence Counter, MicroBeta ² Scintillation & Luminescence Counter	Opera HCS System, Opera LX HCS, Operetta HCS Reader, UltraVIEW [®] VoX	MicroBeta ² Scintillation & Luminescence Counter, TopCount Scintillation & Luminescence Counter, Tri-Carb Liquid Scintillation Analyzer	ViewLux uHTS Imager, EnVision Multilabel Plate Reader with TRF Laser, VICTOR X4 or X5 Multilabel Plate Reader
33, 34, 35, 37	34, 35, 36	31, 37	35	31, 37	35, 36	33, 34

THE LARGEST PORTFOLIO OF GPCR CELL LINES

**Do you need to decrease
your assay development
time and increase the quality
of your results?**

-
- >300 stable GPCR cell lines for binding and functional testing
 - >130 stable GPCR cell lines for calcium testing
 - >150 ready-to-use cAMPZen and AequeZen frozen cells
 - All PerkinElmer products are protected against intellectual property (IP) infringement, unlike those of many competitors
-

PerkinElmer has the largest collection of validated cell lines for GPCR research and offers HTS-friendly cellular testing for agonist and antagonist pharmacology. In addition, the majority of our ValiScreen cell lines have been validated with saturation and competition binding, and approximately 30% of these protocols have been adapted for GTP γ S binding.

- Calcium testing using either aequorin, fluorescent dyes or both
- cAMP testing
- GTP γ S validation
- IP accumulation assays
- AlphaScreen® *SureFire*® ERK validation
- Label-free validation
- Generation of frozen ready-to-use cells

Clone selection and validation

For every PerkinElmer receptor, we offer different stable recombinant clones optimized for specific applications. We have cell lines where clonal selection is carried out with receptor binding. High and low expressors are then tested and selected for optimal binding and functional response according to the receptor-specific transduction pathway. In AequeScreen and PhotoScreen cell lines, clonal selection is carried out with a calcium assay using a luminescent protocol.


Binding conditions are also available for nearly all cell lines. For specific applications on certain receptors, we may select different clones or validate existing clones with alternative protocols. The measure of constitutive activity and its modulation by inverse agonists, for example, often requires high expression levels or adapted protocols.

Cell line validation


Our validated cell lines typically have saturation and competition binding data, along with functional agonist/antagonist dose response curves using reference ligands.

Cell line development

A rigorous cell line development process ensures your cell lines are the highest quality possible.




Radioligand Binding Assay: Adrenergic α_{2B}




Competition binding experiment: binding of [³H]Rauwolschine to membranes expressing the adrenergic α_{2B} receptor was inhibited by increasing concentrations of reference ligands. The separation of the free and bound radioligand was performed by filtration on 96-well filtration plates. Plate was measured on a TopCount. Assay can also be done in homogeneous format using SPA technology.

GTP γ [³⁵S] Binding: Cannabinoid CB₁




GTP γ [³⁵S] binding assay: membranes expressing the Cannabinoid CB₁ receptor were incubated in presence of increasing concentrations of reference ligands and the binding of GTP γ [³⁵S] was measured using the SPA technology. Plate was measured on a TopCount.

LANCE Ultra cAMP (G_i): Cannabinoid CB₁




LANCE Ultra cAMP assay using cAMPZen CHO-CB₁ cells: an antagonist dose-response experiment was performed in 384-well format using 2,500 cells/well. Cell stimulation was performed for 30 min at room temperature, and the agonist (WIN 55,212-2), antagonist and 2.5 μ M Forskolin were added simultaneously. Plate was read on an EnVision.

LANCE Ultra cAMP (G_s): Adrenergic Beta receptor



LANCE Ultra cAMP assay on endogenous β -adrenergic receptors: an agonist dose-response experiment was performed in 384-well format. SK-N-MC cells were used at 2,000 cells/well. Plate was read on an EnVision.

Calcium Flux (Fluo-4): Serotonin 5-HT_{2A}

Ca²⁺ flux assay: cells stably expressing the Serotonin 5-HT_{2A} receptor were incubated in presence of increasing amounts of reference ligand and Ca²⁺ flux were followed using Fluo-4 as tracer. Assay was read on an FDSS 6000 (Hamamatsu Photonics).

A complete offering of cell lines, frozen cells, membrane preparations and complementary ligands.

Fresh Cell Lines


- ValiScreen GPCR and ion channel cell lines
- AequoScreen aequorin recombinant GPCR cell lines, parental cells and plasmids
- PhotoScreen Photina® cell lines

Ready-to-use frozen cells

- AequoZen aequorin frozen cells
- cAMPZen cAMP frozen cells

Complementary Membranes & Ligands

- NEN radioligands
- PerkinElmer membrane preparations

Calcium Flux (Aequorin): Muscarinic M₃

The Muscarinic M₃ AequoScreen cell line agonist (Acetylcholine) and antagonist (N-Me-Scopolamine) dose response curves run on a MicroBeta² LumijET. 384-well format, 6,000 cells/well.

INTRACELLULAR CALCIUM DETECTION MADE EASY

Would easier HTS generic Ca^{2+} measurements make a difference to you?


Measuring intracellular calcium through functional, high throughput screening

of cell lines expressing a GPCR has never been easier. PerkinElmer offers the widest selection of luminescence-based platforms for efficient, cost-effective batch screening. Choose from AequoScreen (aequorin) and PhotoScreen (Photina®) cell lines for validated cells.

We bring you the best performing biological coupling by providing the choice of two different photoproteins to measure calcium signaling, resulting in:


- Large signal-to-noise ratio and low background levels
- Short incubation for tested compounds
- Easily followed reaction kinetics
- Stability of signal unaffected by cell growth

AequoScreen and PhotoScreen Assay Principle



The apophotoprotein requires coelenterazine, to be converted to an active enzyme, aequorin or Photina®, respectively. Upon calcium binding, the photoprotein oxidizes coelenterazine into coelenteramide, producing CO_2 and emitting light. The light emission is measured as luminescence.

CHO Aequorin Vasopressin $\text{V}_{1\text{B}}$ on FLIPR^{TETRA}®



AequoScreen luminescence Ca^{2+} flux assay: cells stably expressing the Vasopressin_{1B} receptor were loaded with coelenterazine h for 4 hours. Ligands were dispensed on cells using the FLIPR^{TETRA} dispensing system and light emission was measured. 384-well adherent assay, 10,000 cells/well, [Arg⁸] Vasopressin $\text{pEC}_{50} = -10.40$.

Shedding new light on calcium measurement

You can be confident in the highest sensitivity of detection regardless of the physiological coupling of the expressed GPCR.

Optimal localization of a Ca^{2+} -activated photoprotein within the mitochondria amplifies the luminescent response triggered by GPCRs coupled to the G_q/PLC pathway. AequoScreen and PhotoScreen deliver high signal-to-background ratio, high throughput (over 200,000 data points per day) and a robust assay (high Z' , low CV values and solvent resistance). They are suitable for both agonist and antagonist screening.


The AequoScreen platform includes:

- Parental cell lines
- Double-transfected GPCR aequorin cell lines
- Frozen ready-to-use cells
- Plasmids

The PhotoScreen platform includes:


- Double-transfected GPCR Photina[®] cell lines

CHO Aequorin H_1 Agonist Assays



Histamine dose-response curve. 1536-well suspension assay, 1,500 cells/well run on a LumiLux[®].

CHO Aequorin H_1 Agonist and Antagonist Assays



384-well agonist and antagonist suspension assays: 5,000 cells/well run on a MicroBeta² LumijET. The dual dispensing feature of the MicroBeta² LumijET microplate reader allows both agonist and antagonist screening in the same well. H_1 AequoScreen cells were dispensed together with either buffer or antagonist (trans-triprolidine, $10 \times \text{pIC}_{50}$) into an OptiPlate™-384. The agonist (histamine, $10 \times \text{pEC}_{50}$) was dispensed and resulting luminescence was recorded (total volume: 70 µL).

CONVENIENT READY-TO-USE FROZEN CELLS


Pay by the assay point

You need to perform aequorin-based calcium measurement or cAMP measurement, but don't have the time or capability for lengthy cell preparation. Now it's easier than ever to have validated, ready-to-use irradiated cells that you just thaw and use. You simply store your cells in a -80 °C freezer or in liquid nitrogen until you're ready to use them.

PerkinElmer AequoZen and cAMPZen ready-to-use frozen cells offer:


- Reliability – cells are validated for assay protocol performance and shipped worldwide
- Flexibility – perform cellular GPCR tests on multiple receptors at a time for screening, lead optimization and profiling; available in aliquots of 2.5 and 10 million cells per vial
- Convenience – you pay for only the cells you need when you need them

Optimized and standardized culture conditions, even for very sensitive GPCRs, eliminate variability in sensitivity and size of response between different batches of cells. Because cells are pre-validated, you avoid testing expensive and time-consuming failed batches. With PerkinElmer frozen cells, you may never have to postpone a scheduled testing campaign again.



Assay window and consistency of EC₅₀ values. The calcium aequorin assay illustrates the correlation between EC₅₀ values of the reference agonist obtained on fresh versus frozen cells with an aequorin readout on 15 receptors. Neurotensin shows the pharmacology comparison values between fresh AequoScreen NTS₁ and equivalent AequoZen frozen cells. Luminescence was measured with the MicroBeta technology.


AequoScreen Assay in CHO-secretin AequoZen Cells



Secretin agonist response in AequoScreen assay. 5,000 cells/well, 384-well format. Luminescence was measured with the MicroBeta technology.

The MicroBeta² LumiJET with multiple injector assemblies and 12-detector assembly is an ideal platform for studying luminescence calcium flux in assay development, high throughput screening and compound screening using aequorin technology.

Stability Testing of Fresh Cells vs. Frozen Cells



Comparison of freshly cultured cells with frozen cell preparations using a cAMP TR-FRET test on MC₄ cells after 15 days and 2 months storage at -80 °C.

MEASURING cAMP? UNMATCHED SENSITIVITY. UNLIMITED POSSIBILITIES.

Choose from several high-performance homogeneous assay technologies to rapidly measure cAMP levels – all designed for optimum detection of $G\alpha_s$ - and $G\alpha_{i/o}$ -coupled GPCRs.


LANCE *Ultra* cAMP: exceptional sensitivity

With 3x greater sensitivity and the widest assay window available, LANCE *Ultra* cAMP offers unsurpassed detection for even the most challenging targets. LANCE *Ultra* cAMP provides more robust data in a shorter run time using 80% fewer cells.

- Optimizes screening of difficult targets and $G\alpha_i$ antagonists
- Delivers robust sensitivity when miniaturized to 1536-well format for μ HTS
- Offers unmatched signal stability – even with overnight incubation
- Uses a simple TR-FRET protocol
- Provides reproducible results with the highest Z' values


LANCE *Ultra* cAMP can be used with PerkinElmer's validated Total GPCR Solution including cAMPZen frozen cell lines.

LANCE *Ultra* cAMP Assay Principle




In the absence of free cAMP (top), maximal TR-FRET signal is achieved. Free cAMP produced by stimulated cells competes with the Eu-cAMP tracer for binding to the ULight™ monoclonal antibody (bottom), causing a decrease in TR-FRET signal proportional to the concentration of cAMP produced.

Sensitivity Unaffected After Overnight Incubation



Unlike other cAMP assay options, the LANCE *Ultra* cAMP assay offers both higher sensitivity and a wider assay window that is not affected by overnight incubation.

Rank Order of Potencies of cAMP Response



Agonist-induced cAMP response in frozen CHO-CXCR3 cells (2,000 cells/well). The observed rank order of potencies and efficacy is consistent with published data (CXCL11 > CXCL10 > CXCL9). Assay was read on an EnVision.

EnVision Multilabel Plate Reader with TRF Laser




The EnVision Multilabel Plate Reader with TRF Laser is ideal for LANCE *Ultra* cAMP assays. This benchtop reader delivers optimized performance for TR-FRET detection and is available with stackers. Optional 20 or 50 microplate magazine stackers make it ideal for use with LANCE *Ultra* cAMP where the same strong TR-FRET signal is observed first plate in to last plate out.

LANCE TR-FRET cAMP: precision detection

Our first-generation LANCE cAMP kit continues to provide outstanding performance in a three-component assay, combining homogeneous TR-FRET technology for maximum excitation/emission discrimination.

- Offers easy assay development and adaptation to automation
- Provides excellent signal-to-background ratios, Z' values and sample throughput for optimum productivity

LANCE cAMP Assay Principle




Light pulse at 340 nm excites the Eu-chelate of the Eu-SA/biotin-cAMP tracer. The energy emitted from the Eu-chelate is transferred to the Alexa Fluor® 647-labeled anti-cAMP antibodies bound to the tracer, generating a TR-FRET signal at 665 nm. Residual energy from the Eu-chelate will produce light at 615 nm. cAMP of a sample competes with the tracer for antibody binding sites and causes signal reduction.

AlphaScreen cAMP assay kit

Exclusive Amplified Luminescent Proximity Homogeneous Assay (Alpha) technology provides an easy-to-automate screening tool for monitoring modulation of adenylyl cyclase through a wide range of $G\alpha_{i/o}$ -coupled and $G\alpha_s$ -coupled GPCRs. This system forms the basis for a high throughput, sensitive, functional cAMP assay often yielding very large signal-to-background ratios and high Z' values.

AlphaScreen cAMP Assay Principle



Detection of cAMP is based on the competition between cAMP produced by cells and a biotinylated cAMP probe recognized by streptavidin-coated Alpha Donor beads and anti-cAMP conjugated Acceptor beads. The beads are brought into proximity and a signal is produced.

**Simple to use cAMP, [¹²⁵I] RIA kit:
scintillation proximity assay**

The homogeneous cAMP [¹²⁵I] RIA kit using SPA beads has been developed for the detection and quantitation of cAMP levels in cell lysates using protein A-coated PVT scintillation beads (RPNQ0019) or protein A-coated PS imaging beads (RPNQ0264). The assay is simple to perform, requiring minimal manipulations, and suitable for automated high throughput screening applications. SPA associated with antibody-bound cAMP can be counted in the presence of unbound radiolabeled cAMP without the need for separation. This homogeneous RIA has a working range of 0.2-12.8 pmol cAMP/well (1.32-84.28 ng/mL) with a sensitivity of detection of 0.1 pmol/well (0.65 ng/mL).

Reliable cAMP, [¹²⁵I] RIA kit: radioimmunoassay

Our easy-to-use cAMP [¹²⁵I] RIA kit measures cAMP in plasma, urine, tissue or cell culture samples. The first and second antibodies are supplied pre-reacted, eliminating an assay step. The assay is based upon competitive ligand binding between iodinated cAMP tracer and cAMP contained in the sample. Depending upon the level of sensitivity required, the assay can be run with acetylation (sensitivity 0.025 pmol/mL) or without.

Rapid FlashPlate [¹²⁵I] assay systems

Adenylyl Cyclase Activation

Solid-phase scintillation proximity FlashPlate microplate technology has long been considered the gold standard for cAMP functional assays in HTS labs. It enables fast and reliable quantification of receptor-mediated adenylyl cyclase activation/inhibition using live cells all in one well. A simple protocol measures a true second messenger response without extraction.

Acetylation/Non-acetylated Standards and Samples

Use FlashPlate cAMP assay kits to measure cAMP in biological fluids, including plasma, urine and tissue samples. Choose acetylated samples when the highest sensitivity is desired.

Functional cAMP Determination

	LANCE <i>Ultra</i>	LANCE	AlphaScreen	cAMP [¹²⁵ I] RIA Kit	cAMP [¹²⁵ I] RIA Kit	FlashPlate [¹²⁵ I]
Platform	TR-FRET	TR-FRET	Luminescence	SPA-RIA	Radioimmunoassay	Scintillation
Homogeneous	√	√	√	√		
Non-rad	√	√	√			
Density Format (well plates)						
96, 384, 1536	√ √	√ √	√ √	√	n/a	√
Throughput (points/day)	>500,000	>500,000	>70,000	>70,000	>500	>20,000
Sample Type						
Whole Cell	√	√	√	√	√	√
Membrane Preparations	√	√	√	√	√	√
Tissue			√	√	√	√
Plasma			√	√	√	√
Urine					√	√
Recommended Microplates	OptiPlate CulturPlate™ ProxiPlate™ 1/2 AreaPlate	OptiPlate CulturPlate ProxiPlate 1/2 AreaPlate	OptiPlate CulturPlate ProxiPlate 1/2 AreaPlate AlphaPlate™	OptiPlate ProxiPlate		FlashPlate
Typical Standard Curve*	IC ₅₀ = 1.4 nM Z' = 0.80	IC ₅₀ = 3 nM Z' = 0.80	IC ₅₀ = 6.0 nM Z' = 0.77	EC ₅₀ = 3 nM	IC ₅₀ = 13.9 nM Z' = 0.98	
>24 Hours Signal Stability	√	√	√	√	√	√

*Comparison of cAMP Assay Technologies for High Throughput Screening. Patricia Kasila and Harry Harney, PerkinElmer, Inc.


THE EASY SOLUTION FOR TOUGH-TO-SCREEN TARGETS

AlphaScreen® SureFire® ERK assay

AlphaScreen® SureFire® is a unique platform for measuring endogenous levels of phosphorylated ERK which provides a powerful tool to measure many GPCR types. It enables the highly sensitive and precise interrogation of various receptors.


- Excellent secondary GPCR technology for challenging targets
- Applicable to endogenous receptors and primary cells
- Alternative for tough-to-screen targets not optimally coupled through the calcium or cAMP pathway

ERK-MAP Kinase and AKT Signaling




ERK is activated by multiple receptor classes, including G_q -, $G_{i/o}$ - and some G_s -coupled receptors, making it an ideal endpoint measurement for GPCR activation.

ERK Phosphorylation in CHO-CCR7 cAMPZen Cells




CHO-CCR7 cAMPZen frozen cells were thawed, seeded at 40,000 cells/well and the next day were stimulated with a dose-response of MIP-3β. A single concentration of ATP, acting on an endogenous P2Y₂ receptor in CHO cells, was used as a control agonist. The cell lysate was used to assess ERK phosphorylation using the AlphaScreen® SureFire® kit, and measured with the EnVision plate reader.

ERK Phosphorylation in CHO-Gal₁ cAMPZen Cells



CHO-Gal₁ cAMPZen frozen cells were thawed, seeded at 40,000 cells/well and the next day were stimulated with a dose-response of galanin. A single concentration of ATP, acting on an endogenous P2Y₂ receptor in CHO cells, was used as a control agonist. The cell lysate was used to assess ERK phosphorylation using the AlphaScreen® SureFire® kit, and measured with the EnVision plate reader.

AlphaScreen® SureFire® ERK Assay Principle



In an AlphaScreen® SureFire® assay, the first antibody is biotinylated and is captured by the streptavidin-coated Alpha Donor bead, which in turn captures the endogenous substrate. The second antibody is captured by the protein A-conjugated Acceptor bead, but only recognizes the phosphorylated form of the substrate of interest. The two beads are only brought into close proximity in the presence of the phosphorylated substrate. This assay setup allows “fishing” of the endogenous substrate in cells, and is extremely sensitive.

Functional ERK Determination

Platform	AlphaScreen® SureFire® Luminescence
Homogeneous	✓
Non-rad	✓
Density Format (well plates) 96, 384, 1536	✓
Throughput (points/day)	>70,000
Sample Type	Most cell types including primary cells & tissue samples
Recommended Microplates	OptiPlates (white, opaque) AlphaPlate ProxiPlate 1/2 AreaPlate

CHOOSE THE RIGHT REPORTER GENE ASSAY FOR YOUR RESEARCH

“Mix-and-measure” luminescence reporter gene assays: **britelite plus, neolite and steadylite plus**

Based on the use of firefly luciferase, these assays provide long-lived luminescent signals that allow efficient processing of functional GPCR targets. Choose either britelite plus, neolite or steadylite plus based on sensitivity and half-life requirements. Each offers:

- Excellent Z' values
- No wash platform
- 2-8 °C storage
- DTT-free – eliminates odor and need for hood work

Strong signaling, remarkably accurate results with **britelite plus**


With an extremely bright signal and a half-life of 30 minutes, the britelite plus homogeneous ultra-high sensitivity reporter gene assay kit detects and clearly highlights a wide range of luciferase content levels in up to 1536-well formats.

neolite – strong and stable

Ideal for medium to high throughput applications that require exceptional sensitivity in batch processing, neolite offers a half-life of at least 2.5 hours with a strong signal intensity. neolite delivers superior reproducibility, as it is less sensitive to mixing conditions than typical luminescent detection systems.

steadylite plus for unequaled signal stability

The steadylite plus system is optimized for use in high-density microplates, reducing the effects of incomplete mixing and capillary action (wicking) often seen in rounded square wells. And the 5-hour half-life will give you the flexibility to run high throughput applications.




neolite with EnSpire multilabel plate reader with ultra-sensitive luminescence shows superior sensitivity at low levels of luciferase enzyme compared to that of a standard luminometer.

Choose the Right Detection Reagent Based on Your Application

	britelite plus	neolite	steadylite plus
Sensitivity	Very High	High	Moderate
Half-life (hours)	0.5	2.5	5
Phenol Red Sensitivity	Moderate	Moderate	Moderate
Plate Format	96, 384, 1536	96, 384, 1536	96, 384, 1536
Low Wicking	√	√	√
Microplates	OptiPlate, CulturPlate, ProxiPlate	OptiPlate, CulturPlate, ProxiPlate	OptiPlate, CulturPlate, ProxiPlate
Applications	<ul style="list-style-type: none"> • Superior sensitivity • Low transfection efficiency • Stem cell transfection • Rapid read time (continuous processing) 	<ul style="list-style-type: none"> • Low transfection efficiency • Stem cell transfection • Increased assay windows • Extended read times 	<ul style="list-style-type: none"> • Strong steady signal • Long extended read times • High throughput screening

Stability Comparison



Stability comparison of lyophilized substrates at 22°C for luciferase detection kits, company P versus britelite plus and steadylite plus, respectively.

ViewLux Ultra HTS Microplate Imager



steadylite plus is able to detect lower levels of luciferase – perfect for more sensitive analyses in 1536-well plates measured on high throughput imagers such as our ViewLux uHTS Imager.


NOW MORE FLEXIBILITY: LABEL-FREE VALIDATION

When you need another option for studying difficult targets or desire more physiologically relevant data, choose from a wide selection of cell lines and frozen cells for label-free technology.

PerkinElmer cell lines and frozen cells offer better flexibility and are fully tested for use with the Corning® Epic® System.


Comparable data can be obtained for both freshly passaged cells and γ -irradiated frozen cells. The pharmacology also compares favorably with labeled technologies.

Comparison of μ -Opioid Receptor Response Profiles: Fresh vs. Frozen Cells



Comparison of response profiles from cells expressing the μ -opioid receptor after stimulation with 15 nM DAMGO (8,000 cells/well, 384-well format; measured on the Epic® reader).

Label-free Assay Performance: Fresh vs. Frozen Cells




The assay performance between freshly passaged cells ($Z' = 0.81$) and γ -irradiated frozen cells ($Z' = 0.83$) also compares favorably.

Benefits of label-free technology

- Ability to study multiple GPCR signaling within a single assay
- Complementary pharmacological profiles compared to labeled technologies (calcium flux, cAMP, GTP γ S, IP, receptor binding)
- Ability to generate more physiologically relevant data
- Ability to study recombinant and endogenous GPCR targets
- No artifacts from working with labels
- Non-invasive
- Pathway independent

Comparison of μ -Opioid Receptor Reference Ligand Efficacy and Potency



Dose response curves obtained with frozen cells (top) EC_{50} for agonist DAMGO and Endomorphin-1 were 2.9 nM and 0.7 nM, respectively. IC_{50} for the antagonist CTOP was 25 nM. Dose response curves obtained with freshly passaged cells (bottom). EC_{50} for agonist DAMGO and Endomorphin-1 were 2.9 nM and 0.9 nM, respectively. IC_{50} for the antagonist CTOP was 36 nM at EC_{90} concentration of the agonist DAMGO (measured with the Epic® reader).


EC₅₀/IC₅₀ Values for the μ -Opioid Receptor are Equivalent for LANCE cAMP and Label-free Technology*

Agonist or Antagonist	EC ₅₀ /IC ₅₀ (nM) Label Free	EC ₅₀ /IC ₅₀ (nM) cAMP
DAMGO	3.2	2.0
Endomorphin	1.0	1.6
CTOP	40.0	40.0

*Freshly passaged cells on the Corning® Epic® System; 2,500 cells/well (cAMP) and 8,000 cells/well (label-free).

Label-free Epic® Assay Principle

The Epic® technology measures the dynamic mass redistribution (DMR) that occurs in response to receptor activation or deactivation. There is an integrated response for every pathway which enables the detection of cellular responses for endogenous as well as over-expressed targets with great sensitivity; the results can be used to complement those obtained using labeled technologies.



Dynamic mass redistribution within a cell causes index of refraction changes resulting in a resonant wavelength shift; the sensing zone of ~150 nm above the sensor is represented by the wave indicated on the left.

MORE REAGENTS TO IMPROVE ASSAY PERFORMANCE

The right mix of GPCR membranes, labeled ligands and consumables to fit every need.

PerkinElmer offers an unsurpassed selection of reagents including radioactive and non-radioactive options to measure GPCR binding. These options can supplement traditional filtration assays to increase the performance and productivity of your assays.

GPCR membrane preparations

Significant Targets: Proven Ligand Receptor Binding

Now you have instant access to the most extensive portfolio of pharmacologically characterized GPCR membrane preparations. Membrane Target Systems® are quality-assured frozen membranes from cells that express recombinant or endogenous receptors.

We submit every batch of receptor to stringent quality control testing that includes saturation radioligand binding to determine receptor density (B_{max}), affinity (K_d) and pharmacological analyses (competition curves). At least 30% of our membrane preparations are also tested for GTP γ S binding.

Ligand receptor binding

DELFLIA Europium-labeled Ligands


The DELFLIA ligand family consists of europium (Eu)-labeled peptides and proteins which are ideal for ligand-receptor binding assay. DELFLIA ligands can be used in several assay formats, including filtration assays using filter plates, as a solid phase assay using streptavidin-coated microplates and biotinylated-WGA or with adherent cells on a cell culture plate.

DELFLIA Eu-labeled ligands are:

- Stable, non-radioactive reagents
- Sensitive – even membranes with low expression levels can be used
- Offered with ready-made protocols for ease of use
- Available for customer labeling

Fluorescent labeled ligands

Fluorescent labeled ligands (CellAura Technologies Ltd.) can also be used as a novel tool with PerkinElmer's recombinant cell line portfolio enabling the determination of ligand binding characteristics using the Opera High Content Screening platform (page 31). Selective membrane binding of the fluorescent ligand can be displaced using unlabeled competitor compounds in addition to being able to study the binding kinetics.



Binding of a fluorescently labeled ligand (β_2 -633AN) to the β_2 adrenergic receptor expressed on HEK293 cells.

NEN radioligands

A pioneer in the development of essential tools for receptor ligand binding, PerkinElmer now offers over 400 NEN ^3H and ^{125}I ligands with new products continuously introduced. Each of our ultra-pure radioligands is fully characterized for pharmacological action and validated in receptor binding assay.


If you need further specialization for a specific application, we will work with you, scientist to scientist, to custom-design a radioligand to meet your needs.

Receptor-ligand binding SPA

The SPA platform technology provides a homogeneous format for studying receptor-ligand binding with either purified receptors or in context, embedded in cellular membranes. SPA beads coated with wheat germ agglutinin (WGA) are used to bind cellular membranes with the appropriate receptor embedded in its normal context or using streptavidin (SA)-coated beads to bind a purified biotinylated receptor. The radiolabeled ligand (^3H -, ^{14}C - or ^{125}I -) bound to the receptor is detected without the need for separation of unbound ligand in suspension.


- Sensitive even membranes with low expression levels can be used
- Extensive portfolio of radiolabeled ligands with matching receptors embedded in cellular membranes
- Available for customer labeling

Saturation Binding Experiment Scintillation Proximity Assay




Saturation binding of [^{125}I]SB 258585 to membranes from 1321N1 cells stably expressing the serotonin 5-HT₆ receptor. Measurement of the bound radioligand was performed with SPA beads.

Saturation Binding Experiment Filtration Method



Saturation binding of [^{125}I]HEAT to membranes from CHO-K1 cells stably expressing the adrenergic α_{1A} receptor. Separation of the bound from the free radioligand was performed on 96-well filtration plate.

Competition Binding Experiment Scintillation Proximity Assay



Competition binding of [^{125}I]SB 258585 to membranes from 1321N1 cells stably expressing the serotonin 5-HT₆ receptor with reference ligands. Measurement of the bound radioligand was performed with SPA beads.

GTP binding assays


Following GPCR activation, GTP replaces GDP on the alpha subunit of the G protein. This causes the alpha unit to disassociate and interact with effectors. PerkinElmer offers radioactive assays to measure GTP exchange.

SPA GTP binding kit

Optimize the conditions for screening novel agonists against a GPCR of interest with our SPA [³⁵S]-GTPγS binding kit (RPNQ0210). This kit is a homogeneous radiometric detection assay based on GDP-GTP exchange on the G protein alpha subunit following GPCR activation by agonists. With the use of this hydrolysis-resistant nucleotide, [³⁵S]-GTPγS, the receptor of interest is kept in an “on” state. Agonist stimulation is measured as an increase in the release of GDP from the G protein and is represented in this assay as an increase in [³⁵S]-GTPγS binding and therefore increased signal.

Exclusive to PerkinElmer, SPA has been widely adopted for many HTS applications due to its extremely high sensitivity and wide dynamic range.

Comparison GTP Binding: SPA vs. Filter Assays



Comparison of A1-agonist stimulated [³⁵S]-GTPγS binding detected by SPA and filter assay. Filter assay with 20 μg protein; SPA with 10 μg protein. Results in dpm (n = 3).


NEN [³⁵S]-GTPγS radionuclide filter binding assay

With PerkinElmer's [³⁵S]-GTPγS and receptor membranes, successful GTP binding assays can be performed in lower throughput with a high degree of sensitivity. This proven performer is efficient and automatable.

FlashPlate technology

PerkinElmer WGA-coated FlashPlate microplates are the foundation of homogeneous assay measurement of the binding of GTP to GPCRs. This assay allows quantitative determination of agonist-induced activation of [³⁵S]-GTPγS binding to isolated GPCR membrane preparations without the need for filtration and wash steps.

Comparison of WGA FlashPlate PLUS to Filtration in Agonist-induced Binding of GTPγS to Human μ Receptor



Human MOR membranes were added to either a 96-well WGA FlashPlate PLUS or a low-protein binding assay plate in the presence of GDP, [³⁵S]-GTPγS and increasing amounts of the μ-opioid-specific agonist DAMGO, and incubated. Contents of low-protein binding assay plate were filtered and counted on the MicroBeta JET.

Selection Criteria for GTP Binding Assays

	Filter Plates	Scintillation Proximity Assay	FlashPlate Microplates
Homogeneous		√	√
Detection Technology	Liquid scintillation	Radiometric scintillation proximity	Radioisotopic scintillation proximity
Density Format			
96-well	√	√	√
384-well		√	√
Throughput (points/day)	>15,000	>70,000	>20,000
Complementary Plate Type	96-well UniFilter® plate	96-well white OptiPlate 384-well white OptiPlate 96-well ProxiPlate 384-well ProxiPlate	WGA FlashPlate microplate
>24 Hours Stability of Signal	√	√	√



PREMIER, POWERFUL HIGH CONTENT SCREENING


Do you want to identify responders on a new level?

High content screening (HCS) and analysis is an invaluable resource for secondary screening and follow-up testing:

- To explore mechanisms of action and assess unwanted side effects of a compound early on
- To follow cAMP or calcium-flux-based primary screens and to provide independent readout to those technologies, including receptor internalization/recycling or arrestin recruitment
- To identify partial agonists and antagonists based on their degree of activating or inhibiting receptor activation; atypical ligands can also be distinguished

This approach is particularly appropriate when measuring the association of the activated GPCR with ancillary signaling molecules such as β -arrestin:

- New methods for mechanisms of actions research
- Functional assays for receptor activation, receptor endocytosis, signal molecule recruitment
- Multiplexed assays with multiparameter readout



In image-based quantification of Endothelin A Receptor (ETAR) Translocation, images were captured on the Opera using a 20x objective. A false color overlay of nuclear (blue) and ETAR (green) fluorescence is shown. Translocation of ETAR becomes apparent as the concentration of the agonist ET-1 is increased. Upon binding of ET-1, the receptor is internalized into pericentriolar recycling endosomes.

UNPARALLELED SCREENING INSIGHT CELL BY CELL

High content analysis can also bring significantly more and different information to calcium-flux-based assays than can whole well readers.

The cell-by-cell-based evaluation provides the



opportunity to identify sub-populations such as high and low responders or responders and non-responders, thus eliminating variations due to transient transfectants and revertants in a population.

Population-based statistics provide readout parameters different from intensity only. Additionally, calcium signaling can be combined with other assay parameters in HCS (multiplexing).

Activation of GPCR calcium release

GPCRs sense molecules outside the cell and activate intracellular signal transduction pathways, which ultimately lead to different cellular responses. Most are linked to calcium signaling. Acapella™ software, which determines sub-populations of cells based on their calcium response, is a particularly useful tool to decrease the number of false positives in drug discovery compared to systems that average the response across the entire population within a well.

CHO-K1 cell lines stably expressing metabotropic glutamate receptors (mGluRs) under the control of an inducible promoter were used in a cell-based calcium flux assay. The intracellular calcium mobilization due to glutamate-induced GPCR activation was detected using the Ca²⁺ responsive dye, Fluo-4. To study calcium flux kinetics, 25 images were captured within 60 seconds and analyzed with the Acapella Kinetic Intensity Analysis script.



(Top) Time course of Fluo-4 mean fluorescence intensity for the whole population (black) and for the active (red) and inactive (blue) sub-populations incubated without MPEP. The intensity was averaged for each time point, over all the cells in each group (231 cells in total). Classification of cells into the two sub-populations was achieved by introducing a peak intensity threshold at 150. (Bottom) An example dose response curve generated from 3 replicate measurements. A significant decrease in the percentage of active cells is observed with increasing concentrations of MPEP.

Opera for High Content Screening



The Opera is a confocal imaging high content screening reader with unmatched resolution and speed. This high-performance system is very flexible to run either fixed cell assays – such as receptor internalization – or live cell kinetic assays such as calcium flux, even in primary screening. High-resolution imaging using high-magnification objectives can detect every detail of a cellular process for every cell in a screen and even follow it over time.

Operetta for High Content Screening



The Operetta is a very easy to use high content screening reader with widefield, brightfield and confocal options. Suited for any size lab, it provides the newcomer to HCS with ready-to-go applications such as receptor internalization/ endocytosis, pit and vesicle formation, translocation and many more. Smaller scale screening campaigns can be run with great level of detail in cellular readout.

UltraVIEW VoX 3D live cell imaging

When a vesicle disappears, where does it go? How does a membrane change across the whole cell? How many bacteria are invading a cell? You face biological mysteries that 2D imaging alone can't unravel. But with the UltraVIEW VoX 3D live cell imaging system, you can find all the answers you're looking for. It's the only 3D spinning disk system that offers acquisition to analysis. Not only will UltraVIEW VoX tell you more about cellular structure and function than other systems, 3D analysis also provides greater context to cellular changes and more accurate results.

The UltraVIEW VoX features leading innovation, such as the advanced Yokogawa® CSU-X1 spinning disk scanner, patented ProSync® technology and award-winning Velocity® software. With this system, you can achieve faster acquisition rates, greater sample protection and higher image data quality than you ever imagined. With the power of 3D live cell imaging, picture all that you can do with your research.

THE RIGHT INSTRUMENTATION. RIGHT WHERE YOU WANT IT

Versatile, scalable PerkinElmer instrumentation solutions perform consistently and reliably no matter what your platform or volume. From stand-alone instrumentation to fully

integrated, automated systems, you can be confident that your GPCR instrument solution has been precisely validated for specific assay technologies. Our application-focused discovery support team means you can be confident our advice will be on target with your needs.

Assay Platform	ViewLux	EnSpire	EnVision	VICTOR X	MicroBeta ²	TopCount	Tri-Carb	WIZARD ²
Alpha Technology		√ (w/Alpha option)	√ (w/Alpha option)					
AequoScreen PhotoScreen			√ (w/ injectors)	√ (w/ dispenser)	√ (MicroBeta ² LumiJET)			
LANCE (TR-FRET)	√		√ (enhanced w/TRF LASER)	√ (X4 & X5 models)				
DELFI (TRF)	√		√ (enhanced w/TRF LASER)	√ (X4 & X5 models)				
Luciferase Assays	√	√	√	√	√	√		
FlashPlate					√	√		
RIA Kits					√	√		√
Radioligands and Radionuclides					√	√	√	√
SPA	√				√	√	√	

EnVision Multilabel Plate Reader

Fast, sensitive and versatile benchtop readers that deliver optimized performance for assay development and high throughput screening. A modular design allows label-specific optical screening.

Choose EnVision to complete your total GPCR Solution for all major assay applications.

- Features bottom-reading, shaking, scanning and kinetics
- Offers temperature control and dispenser options to enable ion channel, calcium flux, dual reporter gene assays and more
- Works with PerkinElmer's AlphaScreen, and TRF (LANCE *Ultra*, DELFIA) for GPCR assays
- Provides scanning and kinetics capabilities and label-specific optical mirror modules and filters to enable GFP assays and dual reporter gene assays using luciferase and beta lactamase
- Performs sensitive kinetic measurements across a wide dynamic range, and high-speed measurements at short-repeat intervals

With easy integration into our automated HCS systems and the ability to accept microplates up to 3456 wells, EnVision's stacker option greatly increases the throughput of kinetic assays: over 100,000 samples per day when 1536-well plates are used.



EnVision Multilabel Plate Reader

EnSpire Multilabel Plate Reader

EnSpire offers affordable detection across multiple platforms including ultra-sensitive luminescence, Alpha Technology and fluorescence/absorbance using quad-monochromator technology. These detection modes enable the measurement of Alpha cAMP, calcium flux and reporter gene assays. With ultra-sensitive luminescence, you can use fewer precious cells to obtain the same robust signal as standard luminescence. Extra sensitivity compensates for low transfection assays, enhancing assay viability. The stronger signal increases your assay window, so there is no need to redesign assays. All this combined can save up to 75% in reduced transfection reagents and substrates.

EnSpire leverages EnVision's world-class performance technology to deliver flexibility and sensitivity, with easy-to-use software, operational temperature control and bottom-reading capability.



EnSpire Multilabel Plate Reader



VICTOR X Multilabel Plate Reader

Multilabel, multitask counters for all light-emitting and light-absorbing detection technologies including:

- Fluorescence (top and bottom)
- Luminescence
- UV absorbance
- Time-resolved fluorometry
- Fluorescence polarization

It is ideal for fast multiwavelength assays in a range of microplate applications, including cell-based assays, toxicology screens and more. With over 10 counting modes, multiple options and a range of detection technologies, including luminescence, the VICTOR X has the right configuration for your assay needs. It can operate as a stand-alone instrument or integrate into robotic systems.

VICTOR X Light Plate Reader

Flexible microplate-based benchtop reader offers considerably more sensitivity for luminescence measurement than standard readers using either aequorin or luciferase assays. It is ideal for a variety of applications including cell-based assays, immunoassays, toxicity screens and gene expression. Use it as a convenient stand-alone system or, with optional accessories, integrate it into robotic systems as needed. It is ideal for PerkinElmer's mix-and-measure luminescence assay technologies.

ViewLux Ultra HTS Microplate Imager

Ultra-high throughput microplate imager for high sensitivity and fast measurement of light from fluorescence polarization (FP), fluorescence intensity, time-resolved fluorescence, luminescence and absorbance assays.

For radiometric detection, the CCD camera in the ViewLux is designed to capture red-shifted light emission from SPA imaging beads in 96-, 384- and 1536-well formats. All HTS screens are rapidly processed with precision in the ViewLux.

Because the instrument reads entire plates in one exposure, throughput is not affected by plate density. Alternate between batch and robot loading according to need – up to 200,000 samples/hr can be read under continuous operation.



ViewLux Ultra HTS Microplate Imager

MicroBeta² and TopCount Liquid Scintillation and Luminescence Counters

The MicroBeta² and TopCount are the most popular scintillation and luminescence counters in the world for low and medium throughput applications. They are ideal for flash luminescence applications such as aequorin/Ca²⁺ measurement, dual-label reporter gene assays and traditional GPCR filtration assays or scintillation proximity assays (SPAs).

If your applications require flash luminescence, choose MicroBeta² LumijET with a choice of single- or dual-channel injectors. Run dual aequorin (agonist and antagonist) screens with no crosstalk between wells for improved data accuracy and assay performance. And with 1-, 2-, 6- or 12-detector configurations, multiple detectors can be used simultaneously to maximize your throughput in formats up to 384-well plates.

All offer linear responses in excess of 20 million CPS and backgrounds below 100 CPS. Temperature control features guarantee the optimum counting conditions for luminescence assays.



TopCount

MicroBeta²



WIZARD² Automatic Gamma Counter

Offering superior counting performance with gamma-emitting samples, the WIZARD² is used for radiometric detection in labs around the world. Featuring Windows[®] XP software, touch-screen operation, easy networking capability and USB ports for simple data transfer, the WIZARD² Automatic Gamma Counter's user-friendly interface system includes a suite of technology advancements for academic, nuclear medicine and pharmaceutical researchers who conduct radiometric immunoassays (RIAs), chromium release studies and positron emission tomography (PET) research among other gamma radionuclide applications.

Don't forget, PerkinElmer has cocktails and vials that were designed for our radiometric instruments. To find the best cocktails and vials for your application, visit www.perkinelmer.com/cocktailselector and use our online cocktails and vials selector tool.



i-Carb Liquid Scintillation Analyzers

The Tri-Carb family of liquid scintillation analyzers (LSA) are benchtop instruments offering the most versatile and sensitive analyzers for detecting small amounts of alpha, beta and gamma radioactivity. Designed for convenient networking, data reduction and data storage, the Tri-Carb family of LSAs offers unmatched flexibility for a variety of applications. A barcode reader option is available for automated sample tracking and the QuantaSmart® instrument and data reduction software includes an enhanced security option for the Tri-Carb 2910TR, 3110TR and 3180TR/SL models.

JANUS Automated Workstation

From compound management and sample preparation to downstream applications for drug discovery, molecular biology and clinical research, we provide innovative solutions for your automated liquid handling needs. Exceptionally easy to use and flexible to grow with your requirements, the JANUS Automated Workstation is easily integrated with our multilabel detection instruments and diverse portfolio of reagents and kits.

Scale from 96- to 384-well plates to improve your assay throughput and eliminate manual errors. And, with multiple pipetting arms and precision dispensing, you won't compromise your assay sensitivity.

Need a more specialized automation solution? Our Integration Solutions Team stands ready to work with you to develop a custom solution that fits your needs.



JANUS Automated Workstation

Opera, Opera LX and Operetta: high resolution, high speed, high content screening

The Opera is the premier confocal microplate imaging reader for simultaneous high-speed and high-resolution screening. Point scan spinning-disc confocal imaging and water immersion lenses generate high resolution, while high speed is attained by high-speed autofocus, simultaneous multicolor image acquisition and multicore parallel data processing. In addition, the Opera is well suited for live cell applications due to very low photo bleaching inherent to the spinning-disc confocal approach. Environmental control and liquid handling can be added.

The Opera LX is an introductory system using the same imaging technology as the Opera. Sequential three-color imaging plus optional non-confocal UV imaging, optional environmental control and liquid handling makes this system almost as flexible and powerful as the full-blown Opera. Both instruments include full Acapella image analysis software including the cellular assay library for online image analysis.

The Operetta is an entry-level instrument ideal as a starter instrument into high content screening. Equipped with a wide variety of fluorescence excitation, emission, brightfield and magnification options, high-speed autofocus plus a confocal unit for background suppression, it is suited to a broad range of image-based assays for receptor activation and signaling pathways. The included

Harmony® image analysis software provides a large number of ready-made applications such as receptor internalization/ endocytosis, pit and vesicle formation or translocation.

Typical applications are secondary screening and follow-up testing with physiologically relevant readout, phenotypic cell assays and RNAi screening at sub-cellular resolution employing live or fixed cell samples.



cell::explorer Automated HCS Platform

To complement the speed and high-resolution imaging capabilities of the Opera and Operetta confocal readers, cell::explorer™ is the only high content screening platform dedicated to complete automation. With robotic plate handling capabilities, the cell::explorer provides true walk-away automation of your cellular screening applications such as primary and secondary screening, RNAi and other compound profiling techniques, without compromising your data quality. With a robust and flexible design, cell::explorer easily supports our portfolio of GPCRs and integrates your choice of microplate reader while dramatically reducing well-to-well variances and improving the handling of your live cells. The combination of our proven expertise and ability to integrate additional instruments makes cell::explorer the premier platform for your cellular screening applications.



Opera High Content Screening System

BETTER MICROPLATES BETTER RESULTS

Precision microplates for every need

Better microplates yield better results. PerkinElmer microplates were designed with both the assay and the instrument in mind for superior screening performance. Miniaturization has never been easier with

our patented 1536-well plate. It's the same size as 96- and 384-well plates, so no instrument or software automation adjustments are necessary. Along with a standard footprint, our 1536-well plate meets the most rigorous performance standards so you can screen with confidence.

No matter the type, size, color or detection technology, plate reader or automated workstation, PerkinElmer offers an extensive selection of high-quality, application-focused microplates to meet the most demanding process needs. Special coating, packaging and barcoding services are also available.



PerkinElmer offers a wide selection of microplates for every application (Top: 96-, 384-, 1536-well OptiPlate; Bottom: 384-well ViewPlate®).

See page 40 for the Microplate Quick Reference Guide by Detection Method or visit www.perkinelmer.com/microplates for our online plate selection tools.

Selection Criteria for Microplates

Product	Application	Format	Assay Platform	Application
OptiPlate	Available in black and white. White polystyrene provides excellent light reflection with highest efficiency and low background for luminescence assays. Black plates can be used with fluorescence detection technologies.	24-, 96-, 384-, 1536-well	Alpha Technologies (AlphaScreen® SureFire® and cAMP), AequoScreen, LANCE, Lites, SPA	cAMP, calcium flux, reporter gene, phospho-ERK, receptor binding
AlphaPlate	Designed specifically for Alpha Technology assays, this light grey plate results in less crosstalk – up to 17 times lower in the 1536-well format	384- and 1536-well, also shallow 384	Alpha Technologies	cAMP, phospho-ERK
ProxiPlate	Available in black and white. Increased signal with shallow well design.	96-, 384-well	Alpha Technologies, LANCE, Lites, SPA	cAMP, reporter gene, phospho-ERK, receptor binding
1/2 AreaPlate	Black or white plates are the same size and well depth as conventional plates with 1/2 the area, 180 µL per well maximum volume	96-well	Alpha Technologies, LANCE	cAMP, phospho-ERK
CulturPlate	Available in black and white. Optimal for cell-based applications providing sterile, tissue-culture treated environment.	24-, 96-, 384-, 1536-well	Alpha Technologies, AequoScreen, LANCE, Lites	cAMP, calcium flux, reporter gene, phospho-ERK
DELFI microplates	Includes clear and yellow plates, which provide low auto-fluorescence	96-Stripwell™ format and 384-well white	DELFI	Receptor binding
FlashPlate	Scintillating coated microplates in either basic format without coating or plus format with coating (streptavidin, wheat germ agglutinin, and others)	96- and 384-well; custom options available	Radioligand binding	Receptor binding, GTPγS
CytoStar-T	Sterile, tissue culture-treated scintillating microplates that are designed for non-invasive live cell assays. The scintillant is integral to the base plate.	96- and 384-well; custom coating options available	Live cell-based assays using SPA technology	SPA-based studies for metabolism; radioligand binding; drug uptake and metabolic turnover; signal transduction
UniFilter filter plate	Designed for ligand binding assays. UniFilter is supplied with either GF/B or GF/C membrane.	96-well	Radioligand binding	Receptor binding, GTPγS
CellCarrier™	High-quality black, treated, clear bottomed plates designed specifically for Opera and Operetta HCS systems	96-, 384- and 1536-well	Imaging	Pathway analysis
PDL/collagen-coated CellCarrier plates	Black clear bottomed plates for cell adhesion and growth	384- and 1536-well	Imaging	Pathway analysis
ViewPlate	Black or white clear bottomed plates for HCS, cell imaging applications and bottom reading	96-, 384- and 1536-well	Imaging, AequoScreen, cell culture	Luminescence/fluorescence calcium flux, imaging, cell culture
PDL/collagen-coated ViewPlates	Black clear bottomed plates for cell adhesion and growth	96-, 384- and 1536-well	Imaging, AequoScreen, cell culture	Luminescence/fluorescence calcium flux, imaging, cell culture
SpectraPlate™	Clear microplates for ELISA and cell culture	96-, 384-well, 384-shallow and 1536-well	Cell culture	

Microplate Quick Reference Guide by Detection Method

Detection Method	Application	Assay Type
Luminescence	Assays	Alpha Technologies (AlphaLISA®, AlphaScreen, AlphaScreen® <i>SureFire</i> ®)
		britelite plus, neolite, steadylite plus
		AequoZen and AequoScreen (Aequorin), PhotoScreen (Photina®)
Fluorescence	Assays	LANCE, LANCE <i>Ultra</i>
		DELFI
	Confocal Imaging	
Radiometric	Filtration Assays	Receptor Binding, Cell Harvesting, ³ H-Thymidine, DNA Binding/Hybrid
	Isotopic Assays	Solid-phase Radiobinding, Scintillation Proximity Assays (SPA)
		Scintillation Proximity Assays (Bead-based SPA)

Visit www.perkinelmer.com/microplates to place an order or find catalog numbers.

PerkinElmer Instrument	Recommended Microplate	Wells/Plate
EnVision with Alpha option EnSpire with Alpha option	AlphaPlate CulturPlate OptiPlate ProxiPlate 1/2 AreaPlate	384, 1536 24, 96, 384, 1536 24, 96, 384, 1536 96, 384 96
EnVision, EnSpire, VICTOR X, ViewLux, TopCount and MicroBeta ²	CulturPlate OptiPlate ViewPlate	24, 96, 384, 1536 24, 96, 384, 1536 96, 384, 1536
EnVision with injectors, EnSpire with injectors, VICTOR X with injectors, MicroBeta ² LumiJET	CulturPlate OptiPlate ViewPlate	24, 96, 384, 1536 24, 96, 384, 1536 96, 384, 1536
EnVision, VICTOR X, ViewLux	CulturPlate OptiPlate ProxiPlate 1/2 AreaPlate	24, 96, 384, 1536 24, 96, 384, 1536 96, 384 96
EnVision, VICTOR X, ViewLux	DELFLA yellow plate DELFLA white plate DELFLA clear plate OptiPlate	96 384 8 x 12 strips 24, 96, 384, 1536
Opera, Operetta	CellCarrier plate ViewPlate	96, 384, 1536 96, 384, 1536
FilterMate™ Harvester, TopCount, MicroBeta ²	Harvest plate OmniFilter plate UniFilter plate	96 96 24, 96
MicroBeta ² , TopCount	Cytostar-T plate FlashPlate ScintiPlate®	96, 384 96, 384 96
MicroBeta ² , TopCount	OptiPlate ProxiPlate	24, 96, 384 96, 384
MicroBeta ²	IsoPlate VisiPlate™	96 24

Validated for

Receptor type	Subtype	Product Type	Host Cell	Binding	GTP-γS	cAMP	Inositol Phosphate	Calcium Luminescence	Calcium Fluorescence	Cat. No.		
GPCRs												
5-Hydroxytryptamine	5-HT _{1A}	AequoScreen*	CHO-K1	X				X		ES-310-A		
		AequoZen*	CHO-K1					X		ES-310-AF		
		ValiScreen†	CHO-K1	X							MCL-508	
				X	X	X					ES-310-C	
		cAMPZen▼	CHO-K1			X					ES-310-CF	
		Membrane preparation	CHO-K1	X								6110501400UA
	HEK293 EBNA		X								RBHS1AM400UA	
	5-HT _{2A}	AequoScreen	CHO-K1	X					X	X	ES-313-A	
		AequoZen	CHO-K1						X		ES-313-AF	
		ValiScreen	CHO-K1	X						X	ES-313-C	
		Membrane preparation	CHO-K1	X							ES-313-M400UA	
	5-HT _{2B}	AequoScreen	CHO-K1	X					X	X	ES-314-A	
		ValiScreen	CHO-K1	X						X	ES-314-C	
		Membrane preparation	CHO-K1	X							ES-314-M400UA	
	5-HT _{2C} non edited isoform	AequoScreen	CHO-K1	X					X	X	ES-318-A	
		AequoZen	CHO-K1						X		ES-318-AF	
		ValiScreen	1321N1	X				X			ES-318-C	
		Membrane preparation	HEK293	X								6110548400UA
			1321N1	X								ES-318-M400UA
	5-HT _{2C} edited isoform	AequoScreen	CHO-K1	X					X	X	ES-315-A	
		AequoZen	CHO-K1						X		ES-315-AF	
		ValiScreen	CHO-K1	X				X		X	ES-315-C	
		Membrane preparation	CHO-K1	X							ES-315-M400UA	
	5-HT _{3A} (ligand-gated ion channel) see Ion Channels											
	5-HT _{5A}	AequoScreen	CHO-K1	X					X	X	ES-401-A	
		AequoZen	CHO-K1						X		ES-401-AF	
		ValiScreen	HEK293	X								MCL-509
			CHO-K1	X								ES-401-C
		Membrane preparation	CHO-K1	X								RBHS5AM400UA
				X								ES-401-M400UA
	5-HT ₆	AequoScreen	CHO-K1	X					X	X	ES-316-A	
		AequoZen	CHO-K1						X		ES-316-AF	
		ValiScreen	CHO-K1/1321N1	X		X					ES-316-C	
cAMPZen		1321N1			X					ES-316-CF		
Membrane preparation		CHO-K1	X								ES-316-M400UA	
	HEK293	X								RBHS6M400UA		
5-HT ₆ (Rat)	Membrane preparation	HEK293	X							RBRS6M400UA		
5-HT ₇	Membrane preparation	CHO-K1	X							6110512400UA		
5-HT ₇ (Rat)	Membrane preparation	HEK293	X							RBRS7M400UA		
Acetylcholine (Muscarinic)	M ₁	AequoScreen	CHO-K1	X				X		ES-210-A		
		AequoZen	CHO-K1					X		ES-210-AF		
		ValiScreen	CHO-K1	X					X	ES-210-C		
		Membrane preparation	CHO-K1	X							RBHM1M400UA	
	M ₂	AequoScreen	CHO-K1	X					X		ES-211-A	
		AequoZen	CHO-K1						X		ES-211-AF	
		ValiScreen	CHO-K1	X	X						ES-211-C	
		Membrane preparation	CHO-K1	X							RBHM2M400UA	

		Validated for									
Receptor type	Subtype	Product Type	Host Cell	Binding	GTP- γ S	cAMP	Inositol Phosphate	Calcium Luminescence	Calcium Fluorescence	Cat. No.	
Adenosine	M ₃	AequoScreen	CHO-K1	X				X		ES-212-A	
		AequoZen	CHO-K1					X		ES-212-AF	
		ValiScreen	CHO-K1	X					X	ES-212-C	
		Membrane preparation	CHO-K1	X						RBHM3M400UA	
	M ₄	AequoScreen	CHO-K1	X				X	X	ES-213-A	
		AequoZen	CHO-K1					X		ES-213-AF	
		ValiScreen	CHO-K1	X	X					ES-213-C	
		Membrane preparation	CHO-K1	X						RBHM4M400UA	
	M ₅	AequoScreen	CHO-K1	X				X	X	ES-214-A	
		AequoZen	CHO-K1	X				X		ES-214-AF	
		ValiScreen	CHO-K1	X			X			ES-214-C	
		Membrane preparation	CHO-K1	X						RBHM5M400UA	
	Adenosine	A ₁	AequoScreen	CHO-K1	X				X	X	ES-010-A
			AequoZen	CHO-K1					X		ES-010-AF
			ValiScreen	CHO-K1	X	X	X				ES-010-C
			cAMPZen	CHO-K1			X				ES-010-CF
			Membrane preparation	CHO-K1	X	X					ES-010-M400UA
		A ₁ (Rat)	Membrane preparation	CHO-K1	X						6110511400UA
				Sf9	X						6110120400UA
		A _{2A}	AequoScreen	HEK293	X				X	X	ES-011-A
AequoZen			HEK293					X		ES-011-AF	
ValiScreen			HEK293	X		X				ES-011-C	
cAMPZen			HEK293			X				ES-011-CF	
Membrane preparation			HEK293	X						RBHA2AM400UA	
A _{2A} (Rat)		ValiScreen	CHO-K1	X						MCL-511	
A _{2B}		AequoScreen	HEK293					X		ES-013-A	
			HEK293	X		X				ES-013-C	
		Membrane preparation	HEK293	X						RBHA2BC	
	HEK293		X			X			ES-013-CF		
	Membrane preparation	HEK293	X						ES-013-M400UA		
A ₃	AequoScreen	CHO-K1	X				X	X	ES-012-A		
	AequoZen	CHO-K1					X		ES-012-AF		
	ValiScreen	CHO-K1	X	X	X				ES-012-C		
	cAMPZen	CHO-K1			X				ES-012-CF		
	PhotoScreen \blacklozenge	CHO-K1					X		AX-001-PCF		
	Membrane preparation	CHO-K1	X	X					ES-012-M400UA		
A ₃ (Rat)	ValiScreen	CHO-K1	X						MCL-512		
	Membrane preparation	HEK293 EBNA	X						RBRA3M400UA		
Adrenoceptors	α_{1A}	AequoScreen	CHO-K1	X				X		ES-036-A	
		AequoZen	CHO-K1					X		ES-036-AF	
		ValiScreen	CHO-K1	X			X			ES-036-C	
		Membrane preparation	CHO-K1	X						ES-036-M400UA	
	α_{1B}	AequoScreen	CHO-K1	X				X	X	ES-037-A	
	α_{1D}	AequoScreen	CHO-K1	X				X		ES-038-A	

* AequoScreen double recombinant photoprotein cell line

 \blacklozenge AequoZen double recombinant photoprotein frozen cell line

† ValiScreen stable recombinant cell line

 \blacktriangledown cAMPZen stable recombinant frozen cell line \blacklozenge PhotoScreen \blacklozenge double recombinant photoprotein cell line

Validated for

Receptor type	Subtype	Product Type	Host Cell	Binding	Validated for						Cat. No.
					GTP- γ S	cAMP	Inositol Phosphate	Calcium Luminescence	Calcium Fluorescence		
	α_{2A}	AequoScreen*	CHO-K1	X				X	X	ES-030-A	
		AequoZen*	CHO-K1					X		ES-030-AF	
		ValiScreen†	MCHO-K1	X	X	X				ES-030-C	
		cAMPZen▼	CHO-K1			X				ES-030-CF	
		Membrane preparation	CHO-K1	X	X					ES-030-M400UA	
			Sf9	X						6110113400UA	
	α_{2B}	AequoScreen	CHO-K1	X				X		ES-031-A	
		AequoZen	CHO-K1					X		ES-031-AF	
		ValiScreen	CHO-K1	X	X					ES-031-C	
		Membrane preparation	CHO-K1	X	X					ES-031-M400UA	
	α_{2C}	AequoScreen	CHO-K1	X				X		ES-032-A	
		ValiScreen	CHO-K1	X		X				ES-032-C	
		cAMPZen	CHO-K1			X				ES-032-CF	
		Membrane preparation	Sf9	X						6110114400UA	
			CHO-K1	X						ES-032-M400UA	
	β_1	AequoScreen	CHO-K1	X				X	X	ES-033-A	
		AequoZen	CHO-K1					X		ES-033-AF	
		ValiScreen	CHO-K1	X		X				ES-033-C	
		cAMPZen	CHO-K1			X				ES-033-CF	
		Membrane preparation	CHO-K1	X						ES-033-M400UA	
		Sf9	X						6110110400UA		
β_2	ValiScreen	CHO-K1	X		X				ES-034-C		
		HEK293 EBNA	X						RBHBE2C		
	cAMPZen	CHO-K1			X				ES-034-CF		
	Membrane preparation	CHO-K1	X						ES-034-M400UA		
		Sf9	X						6110106400UA		
		HEK293	X						RBHBE2M400UA		
β_3	ValiScreen	CHO-K1	X		X				ES-035-C		
	cAMPZen	CHO-K1			X				ES-035-CF		
	Membrane preparation	CHO-K1	X						ES-035-M400UA		
	Non-recombinant, membrane preparation	SK-N-MC	X						RBHBE3M400UA		
Anaphylatoxin	C3a	AequoScreen	CHO-K1	X				X		ES-730-A	
		ValiScreen	CHO-K1	X	X	X				ES-730-C	
		cAMPZen	CHO-K1			X				ES-730-CF	
		Membrane preparation	CHO-K1	X	X					ES-730-M400UA	
	C5a	AequoScreen	CHO-K1	X				X		ES-731-A	
		AequoZen	CHO-K1					X		ES-731-AF	
		ValiScreen	CHO-K1	X	X	X				ES-731-C	
		cAMPZen	CHO-K1			X				ES-731-CF	
	Membrane preparation	CHO-K1	X	X					ES-731-M400UA		
			X						6110526400UA		
Angiotensin	AT ₁	AequoScreen	CHO-K1	X				X	X	ES-072-A	
		AequoZen	CHO-K1					X		ES-072-AF	
		ValiScreen	CHO-K1	X			X			ES-072-C	
		Membrane preparation	Sf9	X						6110121400UA	
		CHO-K1	X							ES-072-M400UA	
	AT ₂	ValiScreen	CHO-K1	X						ES-070-C	
	Membrane preparation	CHO-K1	X						ES-070-M400UA		

		Validated for								
Receptor type	Subtype	Product Type	Host Cell	Binding	GTP- γ S	cAMP	Inositol Phosphate	Calcium Luminescence	Calcium Fluorescence	Cat. No.
	AT ₂ (Mouse)	ValiScreen	CHO-K1	X						ES-071-C
		Membrane preparation	CHO-K1	X						ES-071-M400UA
Apelin	APJ	AequoScreen	CHO-K1					X		ES-460-A
		ValiScreen	CHO-K1	X	X	X				ES-460-C
		cAMPZen	CHO-K1			X				ES-460-CF
		Membrane preparation	CHO-K1	X	X					ES-460-M400UA
Beta Alanine	MGRPRD, TGR7	AequoScreen	CHO-K1				X	X	ES-741-A	
Bombesin	BB ₁	AequoScreen	CHO-K1	X				X	X	ES-581-A
		Membrane preparation	CHO-K1	X						RBHBS1M400UA
	BB ₂	AequoScreen	CHO-K1	X				X	X	ES-582-A
		AequoZen	CHO-K1					X		ES-582-AF
		ValiScreen	CHO-K1	X					X	ES-582-C
		Membrane preparation	HEK293	X						RBHBS2M400UA
		CHO-K1	X						ES-582-M400UA	
	BB ₃	AequoScreen	CHO-K1					X		ES-580-A
Membrane preparation		Balb 3T3	X						RBHBS3M400UA	
Bradykinin	B ₁	AequoScreen	CHO-K1	X				X	X	ES-091-A
		ValiScreen	CHO-K1	X					X	ES-091-C
		Membrane preparation	CHO-K1	X						ES-091-M400UA
	B ₁ (Rat)	AequoScreen	CHO-K1					X		ES-092-A
	B ₂	AequoScreen	CHO-K1	X				X		ES-090-A
		ValiScreen	CHO-K1	X					X	ES-090-C
Membrane preparation		CHO-K1	X						ES-090-M400UA	
Calcitonin	CGRP ₁	AequoScreen	CHO-K1	X				X	X	ES-420-A
		AequoZen	CHO-K1					X		ES-420-AF
		ValiScreen	CHO-K1	X		X				ES-420-C
		cAMPZen	CHO-K1			X				ES-420-CF
		Membrane preparation	CHO-K1	X						ES-420-M400UA
			Sf9	X						6110135400UA
		Non-recombinant, membrane preparation	SK-N-MC	X						RBHGRPM400UA
	AM ₁	Membrane preparation	Sf9	X						6110136400UA
	AM ₂	AequoScreen	CHO-K1	X				X	X	ES-430-A
		ValiScreen	CHO-K1	X		X				ES-430-C
cAMPZen		CHO-K1			X				ES-430-CF	
Membrane preparation		CHO-K1	X						ES-430-M400UA	
Cannabinoid	CB ₁	AequoScreen	CHO-K1					X	X	ES-110-A
		AequoZen	CHO-K1					X		ES-110-AF
		ValiScreen	CHO-K1	X	X	X				ES-110-C
		cAMPZen	CHO-K1			X				ES-110-CF
		Membrane preparation	CHO-K1	X	X					ES-110-M400UA
			Sf9	X						6110129400UA
	HEK293 EBNA	X						RBHCB1M400UA		

* AequoScreen double recombinant photoprotein cell line

◆ AequoZen double recombinant photoprotein frozen cell line

† ValiScreen stable recombinant cell line

▼ cAMPZen stable recombinant frozen cell line

◆ PhotoScreen® double recombinant photoprotein cell line

Validated for

Receptor type	Subtype	Product Type	Host Cell	Validated for						Cat. No.	
				Binding	GTP- γ S	cAMP	Inositol Phosphate	Calcium Luminescence	Calcium Fluorescence		
	CB ₂	ValiScreen†	CHO-K1	X	X	X				ES-111-C	
			HEK293	X						RBXCB2C	
		cAMPZen▼	CHO-K1			X				ES-111-CF	
		Membrane preparation	Sf9	X							6110130400UA
			HEK293 EBNA	X							RBXCB2M400UA
			CHO-K1	X	X					ES-111-M400UA	
Chemerin	GPR1	AequoScreen*	CHO-K1					X		ES-661-A	
		AequoZen♦	CHO-K1					X		ES-661-AF	
Chemokine	CCR1	AequoZen	HEK293	X						RBHCC1C	
		ValiScreen	CHO-K1	X		X				ES-132-C	
cAMPZen		CHO-K1			X				ES-132-CF		
Membrane preparation		HEK293	X							RBHCC1M400UA	
	CCR2b	AequoScreen	CHO-K1	X				X	X	ES-133-A	
		AequoZen	CHO-K1					X		ES-133-AF	
		ValiScreen	CHO-K1	X	X	X				ES-133-C	
		Membrane preparation	CHO-K1	X	X						ES-133-M400UA
			HEK293	X							6110550400UA
	CCR3	AequoScreen	CHO-K1					X		ES-138-A	
		AequoZen	CHO-K1	X				X		ES-138-AF	
		ValiScreen	K562	X		X				ES-138-C	
		Membrane preparation	CHO-K1	X	X					ES-138-M400UA	
	CCR6	AequoScreen	CHO-K1	X				X	X	ES-139-A	
		AequoZen	CHO-K1					X		ES-139-AF	
		ValiScreen	CHO-K1	X	X	X				ES-139-C	
		cAMPZen	CHO-K1			X				ES-139-CF	
		Membrane preparation	CHO-K1	X	X					ES-139-M400UA	
	CCR7	AequoScreen	CHO-K1	X				X	X	ES-140-A	
		AequoZen	CHO-K1					X		ES-140-AF	
		ValiScreen	CHO-K1	X	X	X				ES-140-C	
		cAMPZen	CHO-K1			X				ES-140-CF	
		Membrane preparation	CHO-K1	X	X					ES-140-M400UA	
	CCR8	AequoScreen	CHO-K1					X		ES-136-A	
		AequoZen	CHO-K1					X		ES-136-AF	
		ValiScreen	CHO-K1	X	X	X				ES-136-C	
		cAMPZen	CHO-K1			X				ES-136-CF	
		Membrane preparation	CHO-K1	X	X					ES-136-M400UA	
	CCR9a	AequoScreen	CHO-K1					X		ES-146-A	
		AequoZen	CHO-K1					X		ES-146-AF	
	CCR10	AequoScreen	CHO-K1					X		ES-143-A	
		AequoZen	CHO-K1					X		ES-143-AF	
	CX ₃ CR1	AequoScreen	CHO-K1	X				X	X	ES-137-A	
		AequoZen	CHO-K1					X		ES-137-AF	
		ValiScreen	CHO-K1	X	X	X				ES-137-C	
		cAMPZen	CHO-K1			X				ES-137-CF	
		Membrane preparation	CHO-K1	X	X					ES-137-M400UA	
	CXCR2	AequoScreen	CHO-K1	X				X	X	ES-145-A	
		AequoZen	CHO-K1					X		ES-145-AF	
		ValiScreen	CHO-K1	X	X	X				ES-145-C	
		cAMPZen	CHO-K1			X				ES-145-CF	
		Membrane preparation	CHO-K1	X	X					ES-145-M400UA	
			Sf9	X					6110132400UA		

		Validated for								
Receptor type	Subtype	Product Type	Host Cell	Binding	GTP- γ S	cAMP	Inositol Phosphate	Calcium Luminescence	Calcium Fluorescence	Cat. No.
	CXCR3	AequoScreen	CHO-K1	X				X	X	ES-142-A
		AequoZen	CHO-K1					X		ES-142-AF
		ValiScreen	CHO-K1	X	X	X				ES-142-C
		Membrane preparation	CHO-K1	X	X					ES-142-M400UA
	CXCR6	AequoScreen	CHO-K1					X		ES-720-A
		AequoZen	CHO-K1					X		ES-720-AF
		ValiScreen	CHO-K1		X	X				ES-720-C
		cAMPZen	CHO-K1			X				ES-720-CF
		Membrane preparation	CHO-K1		X					ES-720-M400UA
			HEK293	X						RBHCX6M400UA
	XCR1	AequoScreen	CHO-K1					X		ES-148-A
		AequoZen	CHO-K1					X		ES-148-AF
		ValiScreen	CHO-K1		X	X				ES-148-C
		cAMPZen	CHO-K1			X				ES-148-CF
		Membrane preparation	CHO-K1		X					ES-148-M400UA
Cholecystokinin	CCK ₁	AequoScreen	1321N1	X				X	X	ES-530-A
		AequoZen	1321N1					X		ES-530-AF
		ValiScreen	1321N1	X			X			ES-530-C
		Membrane preparation	1321N1	X						ES-530-M400UA
			Sf9	X						6110125400UA
			NIH-3T3	X						6110508400UA
	CCK ₂	AequoScreen	1321N1	X				X	X	ES-531-A
		AequoZen	1321N1					X		ES-531-AF
		ValiScreen	1321N1	X					X	ES-531-C
		Membrane preparation	1321N1	X						ES-531-M400UA
			HEK293	X						RBHCKBM400UA
Corticotropin-Releasing Factor	CRF ₁	AequoScreen	CHO-K1	X				X	X	ES-152-A
		ValiScreen	CHO-K1	X		X				ES-152-C
		Membrane preparation	HEK293	X						RBHCRF1M400UA
			CHO-K1	X					ES-152-M400UA	
Dopamine	D ₁	ValiScreen	CH4Cl							MCL-505
			L Cells							MCL-506
			CHO-K1	X		X				ES-172-C
		cAMPZen	CHO-K1			X				ES-172-CF
		Membrane preparation	L Cells	X						6110513400UA
			CHO-K1	X						ES-172-M400UA
	D _{2L}	AequoScreen	CHO-K1					X		ES-171-A
		AequoZen	CHO-K1					X		ES-171-AF
		Membrane preparation	Sf9	X						6110137400UA
	D _{2S}	Membrane preparation	CHO-K1	X						RBHD2CM400UA
	D ₃	ValiScreen [†]	CHO-K1	X	X	X				ES-173-C
		Membrane preparation	CHO-K1	X	X					ES-173-M400UA
	D ₃ (Rat)	Membrane preparation	Sf9	X						6110139400UA
D ₄ (Rat)	ValiScreen	CHO-K1	X		X				ES-170-C	

* AequoScreen double recombinant photoprotein cell line

◆ AequoZen double recombinant photoprotein frozen cell line

† ValiScreen stable recombinant cell line

▼ cAMPZen stable recombinant frozen cell line

◆ PhotoScreen® double recombinant photoprotein cell line

Validated for

Receptor type	Subtype	Product Type	Host Cell	Binding	Validated for						Cat. No.
					GTP-γS	cAMP	Inositol Phosphate	Calcium Luminescence	Calcium Fluorescence		
	D _{4.2}	Membrane preparation	Sf9	X						6110112400UA	
			CHO-K1	X						RBHD42M400UA	
	D _{4.4}	ValiScreen†	CHO-K1	X						MCL-507	
		Membrane preparation	Sf9	X						6110122400UA	
	CHO-K1		X						RBHD44M400UA		
	D _{4.7}	Membrane preparation	Sf9	X						6110123400UA	
			CHO-K1	X						RBHD47M400UA	
	Endothelin	ET _A	AequoScreen*	CHO-K1	X				X	X	ES-320-A
			AequoZen♦	CHO-K1					X		ES-320-AF
			ValiScreen	CHO-K1						X	ES-320-C
Membrane preparation			CHO-K1	X						ES-320-M400UA	
ET _B		AequoScreen	CHO-K1	X				X	X	ES-321-A	
		AequoZen	CHO-K1					X		ES-321-AF	
		ValiScreen	CHO-K1	X					X	ES-321-C	
		Membrane preparation	CHO-K1	X						ES-321-M400UA	
Formyl Peptide	FPR2 (FPRL1)	AequoScreen	CHO-K1	X				X	X	ES-610-A	
		AequoZen	CHO-K1					X		ES-610-AF	
		ValiScreen	CHO-K1	X	X	X			X	ES-610-C	
		cAMPZen▼	CHO-K1			X				ES-610-CF	
		Membrane preparation	CHO-K1	X	X					ES-610-M400UA	
Free Fatty Acid	FFA1 (GPR40)	AequoScreen	CHO-K1					X		ES-652-A	
		AequoZen	CHO-K1					X		ES-652-AF	
		ValiScreen	1321N1				X			ES-652-C	
	FFA (GPR120)	AequoScreen	CHO-K1					X		ES-800-A	
		AequoZen	CHO-K1					X		ES-800-AF	
GABA _B	GABA _{B1a}	Membrane preparation	CHO-K1	X						6110545400UA	
			HEK293	X						6110560400UA	
	GABA _{B1b}	Membrane preparation	CHO-K1	X						6110546400UA	
			CHO-K1	X						6110557400UA	
	GABA _{B1a/2}	AequoScreen	CHO-K1	X				X	X	ES-500-A	
		ValiScreen	CHO-K1	X	X	X				ES-500-C	
		cAMPZen	CHO-K1			X				ES-500-CF	
		Membrane preparation	HEK293	X							6110559400UA
CHO-K1			X	X						ES-500-M400UA	
Galanin	GAL ₁	AequoScreen	CHO-K1	X				X		ES-510-A	
		AequoZen	CHO-K1					X		ES-510-AF	
		ValiScreen	CHO-K1	X	X	X				ES-510-C	
		ValiScreen	HEK293 EBNA	X							MCL-513
			HEK293 EBNA	X							RBXGL1C
		cAMPZen	CHO-K1			X				ES-510-CF	
		Membrane preparation	HEK293 EBNA	X							6110537400UA
	CHO-K1		X	X						ES-510-M400UA	
	GAL ₂	AequoScreen	CHO-K1	X				X	X	ES-511-A	
		AequoZen	CHO-K1					X		ES-511-AF	
		ValiScreen	CHO-K1	X	X		X			ES-511-C	
		PhotoScreen♦ Photina double-transfected cell line	CHO-K1					X		AX-008-PCF	
Membrane preparation		CHO-K1	X						ES-511-M400UA		

		Validated for									
Receptor type	Subtype	Product Type	Host Cell	Binding	GTP- γ S	cAMP	Inositol Phosphate	Calcium Luminescence	Calcium Fluorescence	Cat. No.	
Ghrelin	ghrelin	AequoScreen	CHO-K1					X	X	ES-410-A	
		AequoZen	CHO-K1					X		ES-410-AF	
		ValiScreen	CHO-K1	X			X		X	ES-410-C	
		Membrane preparation	HEK293	X							RBHGHS400UA
			CHO-K1	X						ES-410-M400UA	
Glucagon	glucagon	ValiScreen	1321N1	X		X				ES-710-C	
		Membrane preparation	1321N1	X						ES-710-M400UA	
		AequoScreen	CHO-K1					X		ES-710-A	
	GIP	Membrane preparation	HEK293	X						RBHGIP400UA	
	GLP-2	ValiScreen	1321N1			X				ES-711-C	
		cAMPZen	1321N1			X				ES-711-CF	
	Secretin	AequoScreen	CHO-K1					X		ES-712-A	
		AequoZen	CHO-K1					X		ES-712-AF	
Glutamate	mGlu _{5A}	Membrane preparation	CHO-K1	X						ES-555-M400UA	
	mGlu ₇	ValiScreen	CHO-K1			X				ES-553-C	
		cAMPZen	CHO-K1			X				ES-553-CF	
Glycoprotein Hormone	TSH	AequoScreen	CHO-K1					X		ES-790-A	
Gonadotropin-releasing Hormone	GnRH	AequoScreen	CHO-K1	X				X		ES-600-A	
		Membrane preparation	CHO-K1	X						ES-600-M400UA	
GPR	SUCNR1 (Succinate, GPR91)	AequoScreen	CHO-K1					X		ES-744-A	
		AequoZen	CHO-K1					X		ES-744-AF	
	GPR120	AequoScreen	CHO-K1					X		ES-800-A	
		AequoZen	CHO-K1					X		ES-800-AF	
	OXGR1 (GPR99, Citric acid cycle)	AequoScreen	CHO-K1					X		ES-743-A	
		AequoZen	CHO-K1					X		ES-743-AF	
	Histamine	H ₁	AequoScreen	CHO-K1	X				X	X	ES-390-A
			AequoZen	CHO-K1					X		ES-390-AF
ValiScreen			CHO-K1	X					X	ES-390-C	
Membrane preparation			CHO-K1	X						ES-390-M400UA	
H ₂		AequoScreen	CHO-K1	X				X	X	ES-391-A	
		AequoZen	CHO-K1					X		ES-391-AF	
		ValiScreen	CHO-K1	X		X				ES-391-C	
		cAMPZen	CHO-K1			X				ES-391-CF	
		Membrane preparation	HEK293	X						6110565400UA	
		CHO-K1	X							ES-391-M400UA	
H ₃		AequoScreen	CHO-K1	X				X	X	ES-392-A	
		AequoZen	CHO-K1					X		ES-392-AF	
		ValiScreen	CHO-K1	X	X	X				ES-392-C	
		cAMPZen	CHO-K1			X				ES-392-CF	
		Membrane preparation	CHO-K1	X	X					ES-392-M400UA	
H ₄		AequoScreen	CHO-K1	X				X		ES-393-A	
	AequoZen	CHO-K1					X		ES-393-AF		
	ValiScreen	CHO-K1	X	X	X				ES-393-C		
	cAMPZen	CHO-K1			X				ES-393-CF		
	Membrane preparation	CHO-K1	X	X					ES-393-M400UA		

* AequoScreen double recombinant photoprotein cell line

◆ AequoZen double recombinant photoprotein frozen cell line

† ValiScreen stable recombinant cell line

▼ cAMPZen stable recombinant frozen cell line

◆ PhotoScreen® double recombinant photoprotein cell line

Validated for

Receptor type	Subtype	Product Type	Host Cell	Validated for						Cat. No.	
				Binding	GTP- γ S	cAMP	Inositol Phosphate	Calcium Luminescence	Calcium Fluorescence		
KISS1 (Metastin)	KISS1 (GPR54)	AequoScreen*	CHO-K1					X	X	ES-630-A	
		ValiScreen†	MDA-MB-435S	X					X	ES-630-C	
		Membrane preparation	MDA-MB-435S	X							ES-630-M400UA
CHO-K1	X								RBHMTSM400UA		
Leukotriene	BLT ₁ (LTB4R1)	AequoScreen	CHO-K1	X				X	X	ES-340-A	
		ValiScreen	CHO-K1	X	X	X				ES-340-C	
		cAMPZen▼	CHO-K1			X				ES-340-CF	
		Membrane preparation	CHO-K1	X	X					ES-340-M400UA	
	CysLT ₁	AequoScreen	CHO-K1	X				X	X	ES-470-A	
		ValiScreen	CHO-K1	X			X		X	ES-470-C	
		Membrane preparation	CHO-K1	X						ES-470-M400UA	
OXE (HM74-like)	AequoScreen	CHO-K1					X		ES-640-A		
Lysophospholipid	S ₁ P ₂ (EDG5)	AequoScreen	CHO-K1					X		ES-594-A	
		AequoZen♦	CHO-K1					X		ES-594-AF	
	S ₁ P ₄ (EDG6)	AequoScreen	CHO-K1					X	X	ES-592-A	
		AequoZen	CHO-K1					X		ES-592-AF	
	S ₁ P ₃ (EDG8)	AequoScreen	CHO-K1					X	X	ES-593-A	
		AequoZen	CHO-K1					X		ES-593-AF	
Mas related	MrgB3 (Rat)	AequoScreen	CHO-K1					X		ES-745-A	
		AequoZen	CHO-K1					X		ES-745-AF	
	MrgX1	AequoScreen	CHO-K1					X	X	ES-740-A	
	MrgX2	AequoScreen	CHO-K1					X		ES-742-A	
	Melanin-Concentrating Hormone	MCH ₁	AequoScreen	CHO-K1	X				X	X	ES-370-A
AequoZen			CHO-K1					X		ES-370-AF	
ValiScreen			CHO-K1	X	X	X				ES-370-C	
cAMPZen			CHO-K1			X				ES-370-CF	
Membrane preparation			CHO-K1	X	X					ES-370-M400UA	
MCH ₂		AequoScreen	CHO-K1	X				X	X	ES-371-A	
		AequoZen	CHO-K1					X		ES-371-AF	
		ValiScreen	1321N1	X					X	ES-372-C	
		Membrane preparation	1321N1	X						ES-372-M400UA	
Melanocortin		MC ₁	ValiScreen	CHO-K1	X		X				ES-195-C
	cAMPZen		CHO-K1			X				ES-195-CF	
	Membrane preparation		CHO-K1	X						ES-195-M400UA	
	MC ₃	ValiScreen	CHO-K1	X		X				ES-193-C	
		cAMPZen	CHO-K1			X				ES-193-CF	
		Membrane preparation	HEK293	X							RBXMC3M400UA
	CHO-K1		X							ES-193-M400UA	
	MC ₃ (Mouse)	ValiScreen	CHO-K1	X		X				ES-190-C	
		cAMPZen	CHO-K1			X				ES-190-CF	
		Membrane preparation	CHO-K1	X						ES-190-M400UA	
	MC ₄	AequoScreen	CHO-K1					X	X	ES-191-A	
		AequoZen	CHO-K1					X		ES-191-AF	
		ValiScreen	CHO-K1	X		X				ES-191-C	
		cAMPZen	CHO-K1			X				ES-191-CF	
		Membrane preparation	HEK293	X							RBHMC4M400UA
			CHO-K1	X							ES-191-M400UA

Validated for										
Receptor type	Subtype	Product Type	Host Cell	Binding	GTP- γ S	cAMP	Inositol Phosphate	Calcium Luminescence	Calcium Fluorescence	Cat. No.
	MC ₅	ValiScreen	CHO-K1	X		X				ES-194-C
			HEK293 EBNA	X						RBXMC5C
		cAMPZen	CHO-K1			X				ES-194-CF
		Membrane preparation	HEK293 EBNA	X						RBXMC5M400UA
				CHO-K1	X					ES-194-M400UA
	MC ₅ (Mouse)	ValiScreen	CHO-K1	X		X				ES-192-C
		cAMPZen	CHO-K1			X				ES-192-CF
		Membrane preparation	CHO-K1	X						ES-192-M400UA
Melatonin	MT ₁	AequoScreen	CHO-K1	X				X		ES-620-A
		AequoZen	CHO-K1					X		ES-620-AF
		ValiScreen	CHO-K1	X	X					ES-620-C
		Membrane preparation	CHO-K1	X	X					ES-620-M400UA
	MT ₂	AequoScreen	CHO-K1	X				X		ES-621-A
		AequoZen	CHO-K1	X				X		ES-621-AF
		ValiScreen	CHO-K1	X	X					ES-621-C
		Membrane preparation	CHO-K1	X	X					ES-621-M400UA
Motilin	motilin	AequoScreen	CHO-K1	X				X	X	ES-380-A
		ValiScreen	HEK293	X						RBHMOTC
			CHO-K1	X			X		X	ES-380-C
		Membrane preparation	HEK293	X						RBHMOTM400UA
			CHO-K1	X						ES-380-M400UA
Muscarinic (see Acetylcholine)										
Neurokinin (see Tachykinin)										
Neuromedin U	NMU1	AequoScreen	CHO-K1	X				X	X	ES-450-A
		ValiScreen	CHO-K1	X			X		X	ES-450-C
		Membrane preparation	HEK293	X						RBHNU1M400UA
			CHO-K1	X						ES-450-M400UA
	NMU2	AequoScreen	CHO-K1	X				X	X	ES-451-A
		ValiScreen	CHO-K1	X			X		X	ES-451-C
		Membrane preparation	HEK293	X						RBHNU2M400UA
			CHO-K1	X						ES-451-M400UA
Neuropeptide S	NPS	AequoScreen	CHO-K1	X				X	X	ES-770-A
		ValiScreen	CHO-K1	X			X		X	ES-770-C
Neuropeptide FF	NPFF1	ValiScreen transfected cell line	CHO-K1	X		X				ES-491-C
		cAMPZen	CHO-K1			X				ES-491-CF
		Membrane preparation	CHO-DHRF(-)	X						RBHNF1M400UA
			CHO-K1	X						ES-491-M400UA
	NPFF2	AequoScreen	CHO-K1	X				X	X	ES-490-A
		ValiScreen	CHO-K1	X	X	X				ES-490-C
		cAMPZen	CHO-K1			X				ES-490-CF
		Membrane preparation	CHO-K1	X						RBHNF2M400UA
Neuropeptide W/B	NPBW2 (GPR8)	Membrane preparation	HEK293	X					RBHGP8M400UA	

* AequoScreen double recombinant photoprotein cell line
 ♦ AequoZen double recombinant photoprotein frozen cell line
 † ValiScreen stable recombinant cell line

▼ cAMPZen stable recombinant frozen cell line
 ♦ PhotoScreen® double recombinant photoprotein cell line

Validated for

Receptor type	Subtype	Product Type	Host Cell	Binding	Validated for						Cat. No.
					GTP-γS	cAMP	Inositol Phosphate	Calcium Luminescence	Calcium Fluorescence		
Neuropeptide Y	Y ₁	AequoScreen*	CHO-K1	X				X		ES-351-A	
		ValiScreen†	CHO-K1	X		X				ES-351-C	
		cAMPZen▼	CHO-K1			X				ES-351-CF	
		Membrane preparation	Sf9	X						6110133400UA	
		Non-recombinant, membrane preparation	SK-N-MC	X						RBHNP1M400UA	
	Y ₂	AequoScreen	CHO-K1	X				X		ES-352-A	
		AequoZen*	CHO-K1					X		ES-352-AF	
		ValiScreen	CHO-K1	X		X				ES-352-C	
		cAMPZen	CHO-K1			X				ES-352-CF	
		Non-recombinant, membrane preparation	KAN-TS	X						RBHNP2M400UA	
Neurotensin	NTS ₁	AequoScreen	CHO-K1	X				X	X	ES-690-A	
		AequoZen	CHO-K1					X		ES-690-AF	
		ValiScreen	CHO-K1	X					X	ES-690-C	
			HEK293	X						RBXNT1C	
		Membrane preparation	CHO-K1	X						6110518400UA	
			HEK293	X						RBXNT1M400UA	
	CHO-K1		X						ES-690-M400UA		
	NTS ₂	ValiScreen	1321N1	X						ES-691-C	
		Membrane preparation	HEK293	X						6110566400UA	
			1321N1	X						ES-691-M400UA	
Nicotinic	GPR109A (HM74A)	AequoScreen	CHO-K1	X				X	X	ES-760-A	
		ValiScreen	CHO-K1	X	X	X				ES-760-C	
		cAMPZen	CHO-K1			X				ES-760-CF	
		Membrane preparation	CHO-K1	X	X					ES-760-M400UA	
Opioid	delta	Membrane preparation	CHO-K1	X						RBHODM400UA	
			HEK293	X						6110549400UA	
	delta (Mouse)	Membrane preparation	Sf9	X						6110115400UA	
	kappa	AequoScreen	CHO-K1	X				X	X	ES-541-A	
		ValiScreen	CHO-K1	X	X	X				ES-541-C	
		Membrane preparation	HEK293	X						6110558400UA	
			CHO-K1	X	X					ES-541-M400UA	
	mu	AequoScreen	CHO-K1	X				X	X	ES-542-A	
		AequoZen	CHO-K1					X		ES-542-AF	
		ValiScreen	CHO-K1	X	X	X				ES-542-C	
		cAMPZen	CHO-K1			X				ES-542-CF	
		Membrane preparation	CHO-K1	X	X					ES-542-M400UA	
	NOP (ORL1)	AequoScreen	CHO-K1	X				X	X	ES-230-A	
		ValiScreen	CHO-K1	X	X	X				ES-230-C	
			HEK293	X						RBHORLC	
cAMPZen		CHO-K1			X				ES-230-CF		
Membrane preparation		CHO-K1	X						6110540400UA		
		HEK293	X						RBHORLM400UA		
Orexin	OX ₁	AequoScreen	CHO-K1	X				X	X	ES-330-A	
		AequoZen	CHO-K1					X		ES-330-AF	
		ValiScreen	CHO-K1	X					X	ES-330-C	
		Membrane preparation	CHO-K1	X						ES-330-M400UA	

		Validated for								
Receptor type	Subtype	Product Type	Host Cell	Binding	GTP- γ S	cAMP	Inositol Phosphate	Calcium Luminescence	Calcium Fluorescence	Cat. No.
	OX ₂	AequoScreen	CHO-K1	X				X	X	ES-331-A
		AequoZen	CHO-K1					X		ES-331-AF
		ValiScreen	CHO-K1	X					X	ES-331-C
Peptide P518	QRFP (OX ₂ -like, GPR103)	AequoScreen	CHO-K1	X				X	X	ES-670-A
		ValiScreen	CHO-K1	X					X	ES-670-C
		Membrane preparation	CHO-K1	X						ES-670-M400UA
			HEK293	X						RBH103M400UA
Prokineticin	PKR ₁	AequoScreen	CHO-K1	X				X	X	ES-750-A
		AequoZen	CHO-K1					X		ES-750-AF
		Membrane preparation	CHO-K1	X						ES-750-M400UA
	PKR ₂	AequoScreen	CHO-K1	X				X	X	ES-751-A
		AequoZen	CHO-K1					X		ES-751-AF
		Membrane preparation	CHO-K1	X						ES-751-M400UA
			HEK293	X						RBHPK2M400UA
		ValiScreen	CHO-K1			X				ES-751-C
		cAMPZen	CHO-K1			X				ES-751-CF
Prolactin Releasing Peptide	PRRP	AequoScreen	CHO-K1	X				X	X	ES-480-A
		ValiScreen	CHO-K1	X			X			ES-480-C
		Membrane preparation	CHO-K1	X						ES-480-M400UA
Prostanoid	DP ₁	AequoScreen	1321N1	X				X	X	ES-560-A
		ValiScreen	1321N1	X		X				ES-560-C
	DP ₂ (CRTH2)	AequoScreen	CHO-K1	X				X	X	ES-561-A
		ValiScreen	CHO-K1	X	X	X				ES-561-C
		cAMPZen	CHO-K1			X				ES-561-CF
		Membrane preparation	HEK293	X						RBH2M400UA
	CHO-K1		X	X					ES-561-M400UA	
	EP ₂	ValiScreen	HEK293	X		X				ES-562-C
		cAMPZen	HEK293			X				ES-562-CF
		Membrane preparation	HEK293	X						ES-562-M400UA
	EP ₄	AequoScreen	HEK293	X				X	X	ES-563-A
		ValiScreen	HEK293	X		X				ES-563-C
	FP	ValiScreen	1321N1	X					X	ES-564-C
	TP	ValiScreen	HEK293 EBNA	X						RBHTPC
		Membrane preparation	HEK293 EBNA	X						RBHTPM400UA
Protease-activated	PAR2	AequoScreen	CHO-K1					X	X	ES-780-A
	PAR4	AequoScreen	CHO-K1					X		ES-782-A
		AequoZen	CHO-K1					X		ES-782-AF
Purinergic P2Y	P2Y ₁	ValiScreen	1321N1	X						RBHP1C
	P2Y ₂	ValiScreen	1321N1				X			ES-290-C
	P2Y ₄	ValiScreen	1321N1	X						RBHP4C
	P2Y ₆	ValiScreen	1321N1				X			ES-292-C
	P2Y ₁₁	AequoScreen	1321N1					X		ES-293-A
		AequoZen	1321N1					X		ES-293-AF
		ValiScreen	1321N1	X						RBHP11C
	P2Y ₁₂	ValiScreen	CHO-K1	X		X				RBHP12C

* AequoScreen double recombinant photoprotein cell line

◆ AequoZen double recombinant photoprotein frozen cell line

† ValiScreen stable recombinant cell line

▼ cAMPZen stable recombinant frozen cell line

◆ PhotoScreen® double recombinant photoprotein cell line

Validated for

Receptor type	Subtype	Product Type	Host Cell	Validated for						Cat. No.	
				Binding	GTP- γ S	cAMP	Inositol Phosphate	Calcium Luminescence	Calcium Fluorescence		
Relaxin	RXFP3	AequoScreen*	CHO-K1					X		ES-656-A	
		ValiScreen†	CHO-K1		X	X				ES-656-C	
		cAMPZen▼	CHO-K1			X				ES-656-CF	
		Membrane preparation	CHO-K1		X					ES-656-M400UA	
	RXFP4	ValiScreen	CHO-K1		X	X				ES-659-C	
		Membrane preparation	CHO-K1		X					ES-659-M400UA	
AequoScreen		CHO-K1					X		ES-659-A		
Serotonin (see 5-hydroxytryptamine)											
Somatostatin	sst ₁	ValiScreen	CHO-K1	X	X					ES-520-C	
		Membrane preparation	CHO-K1	X	X					ES-520-M400UA	
	sst _{2a}	AequoScreen	CHO-K1	X				X	X	ES-521-A	
		AequoZen♦	CHO-K1					X		ES-521-AF	
		ValiScreen	CHO-K1	X	X	X				ES-521-C	
		cAMPZen	CHO-K1			X				ES-521-CF	
		Membrane preparation	CHO-K1	X	X					ES-521-M400UA	
		sst ₃	AequoScreen	CHO-K1	X				X	X	ES-523-A
	ValiScreen		CHO-K1	X	X	X				ES-523-C	
	cAMPZen		CHO-K1			X				ES-523-CF	
	Membrane preparation		CHO-K1	X	X					ES-523-M400UA	
	sst ₄	AequoScreen	CHO-K1	X				X	X	ES-524-A	
		AequoZen	CHO-K1					X		ES-524-AF	
		ValiScreen	HEK293 EBNA	X						RBHST4C	
			CHO-K1	X	X	X				ES-524-C	
		cAMPZen	CHO-K1			X				ES-524-CF	
		Membrane preparation	HEK293 EBNA	X						RBHST4M400UA	
	sst ₅			CHO-K1	X	X				ES-524-M400UA	
		AequoScreen	CHO-K1	X				X	X	ES-522-A	
			CHO-K1					X		ES-522-AF	
		ValiScreen	HEK293 EBNA	X						RBHST5C	
			CHO-K1	X	X	X				ES-522-C	
		cAMPZen	CHO-K1			X			ES-522-CF		
		Membrane preparation	HEK293 EBNA	X						RBHST5M400UA	
			CHO-K1	X	X					ES-522-M400UA	
		Succinate (see GPR91)									
		Tachykinin	NK ₁	Non-recombinant, membrane preparation	UC11	X					
NK ₂	AequoScreen			CHO-K1	X				X	X	ES-251-A
	ValiScreen		CHO-K1	X					X	ES-251-C	
	Membrane preparation		CHO-K1	X						ES-251-M400UA	
NK ₃	AequoScreen		CHO-K1	X				X	X	ES-252-A	
	ValiScreen		CHO-K1	X			X			ES-252-C	
	Membrane preparation		CHO-K1	X						ES-252-M400UA	
Thyrotropin Releasing Hormone	TRH ₁		AequoScreen	CHO-K1	X				X	X	ES-700-A
			AequoZen	CHO-K1					X		ES-700-AF
		ValiScreen	CHO-K1	X					X	ES-700-C	
		Membrane preparation	CHO-K1	X						ES-700-M400UA	
Urotensin	UT (Uro II, GPR14)	AequoScreen	CHO-K1	X				X	X	ES-440-A	
		ValiScreen	CHO-K1	X					X	ES-440-C	
		PhotoScreen♦ Photina double-transfected cell line	CHO-K1					X		AX-007-PCF	
		Membrane preparation	CHO-K1	X						ES-440-M400UA	

				Validated for						
Receptor type	Subtype	Product Type	Host Cell	Binding	GTP- γ S	cAMP	Inositol Phosphate	Calcium Luminescence	Calcium Fluorescence	Cat. No.
Vasoactive Intestinal Peptide	UT (Uro II, GPR14) (Rat)	AequoScreen	CHO-K1	X				X		ES-441-A
		AequoZen	CHO-K1					X		ES-441-AF
	UT (Uro II, GPR14) (Mouse)	Membrane preparation	CHO-K1	X						RBMUR2M400UA
	PAC ₁ (PACAP)	AequoScreen	CHO-K1	X				X	X	ES-272-A
		VPAC ₁ (VIP)	AequoScreen	CHO-K1	X				X	X
VPAC ₂	Non-recombinant, membrane preparation	HT29	X							RBHV1PM400UA
	AequoScreen	CHO-K1					X	X	ES-274-A	
Vasopressin	V _{1A}	AequoScreen	CHO-K1	X				X	X	ES-361-A
		AequoZen	CHO-K1					X		ES-361-AF
		ValiScreen	131N1	X					X	ES-361-C
		Membrane preparation	HEK293	X						
	131N1		X							ES-361-M400UA
	V _{1B}	AequoScreen	CHO-K1	X				X	X	ES-362-A
		AequoZen	CHO-K1					X		ES-362-AF
		ValiScreen	131N1	X					X	ES-362-C
			HEK293	X						RBHV1BC
		Membrane preparation	CHO-K1	X						6110543400UA
			HEK293	X						RBHV1BM400UA
	131N1		X						ES-362-M400UA	
	V _{1B} (Rat)	ValiScreen	CHO-K1	X					X	ES-360-C
		Membrane preparation	CHO-K1	X						ES-360-M400UA
	V ₂	AequoScreen	CHO-K1	X				X	X	ES-363-A
		ValiScreen	131N1	X		X				ES-363-C
		cAMPZen	131N1			X				ES-363-CF
		Membrane preparation	MDA-MB-435	X						6110541400UA
			131N1	X						ES-363-M400UA
Tyrosine Kinase Receptors										
Epidermal Growth Factor	EGF	Non-recombinant, membrane preparation	A431	X						RBHEGFM400UA
Transporters										
Dopamine	n.a.	Membrane preparation	CHO-K1	X						RBHDATM400UA
Norepinephrine	n.a.	Membrane preparation	MDCK	X						RBHNETM400UA
Serotonin	n.a.	Membrane preparation	HEK293	X						RBHSTM400UA
Ion Channels										
5-Hydroxytryptamine	5-HT _{3A}	AequoScreen	HEK293					X		ES-402-A
		AequoZen	HEK293					X		ES-402-AF
	5-HT ₃	Membrane preparation	HEK293	X						RBHS3M400UA
Calcium Release Activated	CRAC (Stim1 + Orai1)	PhotoScreen validated in e-phys	HEK293				X		AX-013-PCL	
hERG	K ⁺ channel	Membrane preparation	HEK293	X						RBHERGM400UA
Nicotinic	n.a.	Non-recombinant, membrane preparation	IRM32	X						RBHNICM400UA
Purinergic	P2RX ₂ / P2RX ₃ (Human/Rat)	PhotoScreen validated in e-phys	CHO-K1					X		AX-014-PCL
	P2RX ₁	PhotoScreen validated in e-phys	CHO-K1					X		AX-009-PCF
	P2RX ₄	PhotoScreen validated in e-phys	HEK293					X		AX-015-PCF

* AequoScreen double recombinant photoprotein cell line

◆ AequoZen double recombinant photoprotein frozen cell line

† ValiScreen stable recombinant cell line

▼ cAMPZen stable recombinant frozen cell line

◆ PhotoScreen® double recombinant photoprotein cell line

Reference Guide and Ordering Information

Receptor type	Subtype	Product Type	Validated for						Cat. No.	
			Binding	GTPγS	cAMP	Inositol Phosphate	Calcium Luminescence	Calcium Fluorescence		
Transient receptor potential	TRPA1	PhotoScreen♦ validated in e-phys	HEK293					X		AX-004-PCL
	TRPC3	ValiScreen† validated in e-phys and membrane potential assay	HEK293							AX-011-C
	TRPC6	ValiScreen validated in e-phys and membrane potential assay	HEK293							AX-012-C
Voltage gated potassium channel	Kv1.3 (KCNA3)	ValiScreen validated in e-phys and membrane potential assay	CHO-DUKX							AX-010-C
Wild-type membranes										
1321N1 cells		Membrane preparation	1321N1							RBH1321M010MG
A9 cells	Mouse	Membrane preparation	A9							RBMA9M010MG
CHO-K1 cells	Hamster	Membrane preparation	CHO-K1							RBCCHOM010MG
HEK293 cells		Membrane preparation	HEK293							RBHHEKM010MG
K562 cells		Membrane preparation	K562							RBHK56M010MG
THP-1 cells		Membrane preparation	THP-1							RBHTHPM010MG
U373 cells		Membrane preparation	U373							RBHU37M010MG
Parental cell lines										
		AequoScreen*	1321N1					X		ES-000-A27
		AequoZen♦	1321N1					X		ES-000-A27F
		AequoScreen	1321N1					X		ES-000-A28
		AequoZen	1321N1					X		ES-000-A28F
		AequoScreen (+Gα _{1β} , medium level)	CHO-K1					X		ES-000-A2
	Hamster	AequoZen s (+Gα _{1β} , medium level)	CHO-K1					X		ES-000-A2F
	Hamster	AequoScreen	CHO-K1						X	ES-000-A12
	Hamster	AequoZen	CHO-K1					X		ES-000-A12F
	Hamster	AequoScreen	CHO-K1					X		ES-000-A21
	Hamster	AequoZen	CHO-K1					X		ES-000-A21F
	Hamster	AequoScreen	CHO-K1					X		ES-000-A24
	Hamster	AequoZen	CHO-K1					X		ES-000-A24F
		AequoScreen	HEK293					X		ES-000-A30
		AequoZen	HEK293					X		ES-000-A30F
		AequoScreen	HEK293					X		ES-000-A26
		AequoZen	HEK293					X		ES-000-A26F
Aequorin Plasmids										
		AequoScreen						X		ES-002-AC
		AequoScreen						X		ES-003-AC
Starter Kits										
Muscarinic and Histamine	M ₃ and H ₁	AequoZen Starter Kit	CHO-K1					X		ES-001-AF
Muscarinic and Histamine	M ₃ and H ₁	AequoScreen Starter Kit	CHO-K1					X		ES-001-A

* AequoScreen double recombinant photoprotein cell line

♦ AequoZen double recombinant photoprotein frozen cell line

† ValiScreen stable recombinant cell line

▼ cAMPZen stable recombinant frozen cell line

◆ PhotoScreen® double recombinant photoprotein cell line

PerkinElmer – Your Global Discovery Partner

PerkinElmer gives you easy access to a worldwide network of local support offices, staffed by expert professionals who speak your language. You'll get expert assistance for all your questions – no matter where you are located. Contact us at **800-762-4000**, or visit www.perkinelmer.com/ContactUs.

Keeping up to date on the ever changing range of GPCR research solutions is easy. Simply visit www.perkinelmer.com/GPCRcomplete.

PerkinElmer, Inc.
940 Winter Street
Waltham, MA 02451 USA
P: (800) 762-4000 or
(+1) 203-925-4602
www.perkinelmer.com



For a complete listing of our global offices, visit www.perkinelmer.com/ContactUs

Copyright ©2009-2010, PerkinElmer, Inc. All rights reserved. PerkinElmer® is a registered trademark of PerkinElmer, Inc. *SureFire* is a registered trademark of TGR BioSciences Pty Ltd. All other trademarks are the property of their respective owners.