



PRODUCT DATA SHEET

CYTO-ID[®] Green long-term cell tracer kit

ENZ-51036

Live cell fluorescent labeling over extended time periods with no apparent toxic effects

Product Number/Sizes

ENZ-51036-K025

1 Kit

- Allows dual labeling with a variety of CELLESTIAL[®] fluorescent probes
- Minimal transfer of fluorescence from dye-labeled to unlabeled cells
- Suitable for long-term cell viability, cytotoxicity, cell adhesion, cell migration and cell-cell fusion assays

CYTO-ID[®] Green long-term cell tracer kit uses proprietary noncovalent cell labeling technology to stably incorporate a green fluorescent dye containing hydrophobic aliphatic chains into the cell membrane's lipid bilayer. The dye may be loaded into cells by following the included protocol. The labeling buffer is isotonic for mammalian cells and contains no detergents or organic solvents. The appearance of labeled cells may vary depending upon the cell type from uniformly bright to punctuate. This difference is thought to relate to the extent of membrane internalization occurring after cell labeling. The CYTO-ID[®] Green tracer dye fluorescence is independent of pH within normally encountered physiologic ranges and fluorescence intensity per cell is typically unaffected by the ultimate pattern of dye distribution. The CYTO-ID[®] Green tracer dye is not toxic to cells, as determined using the benchmark MTT cell viability assay. The dye is well retained by cells for up to 96 hours after loading, and is passed to daughter cells upon mitosis. Since the dye does not covalently modify proteins within the cells, normal physiological responses are better preserved than with molecular probes based upon thiol-reactive chloromethyl-based or amine-reactive succinimidyl ester-based fluorescent dyes. Dual labeling is also possible using a variety of available CELLESTIAL[®] dyes. Labeled cells can be visualized by epifluorescence or confocal fluorescence microscopy. Additionally, dye-labeled and unlabeled cell populations can be analyzed by flow cytometry. No transfer of fluorescence to adjacent cells was observed after a prolonged 96 hour incubation period. This is in stark contrast to Calcein AM and BCECF AM, which are only retained within viable cells for a few hours at physiological temperatures. The kit is suitable for a variety of applications including long term cell viability, cytotoxicity, cell adhesion, cell migration and cell-cell fusion studies.

Product Specifications

APPLICATIONS:	Flow Cytometry
APPLICATION NOTES:	CYTO-ID [®] Green long-term cell tracer kit is suitable for a variety of applications including long term cell viability, cytotoxicity, cell adhesion, cell migration and cell-cell fusion studies.
QUALITY CONTROL:	A sample from each lot of CYTO-ID [®] Green long-term cell tracer kit is used to stain Jurkat cells and analyzed by flow cytometry, using the procedures described in the user manual. Mean fluorescence of stained to unstained cells is greater than 5.
QUANTITY:	25 assays
HANDLING:	Protect from light. Avoid freeze/thaw cycles.
SHIPPING:	Shipped on Dry Ice
SHORT TERM STORAGE:	-20°C
LONG TERM STORAGE:	-20°C
KIT/SET CONTAINS:	CYTO-ID [®] Green tracer dye, 50µl 4X Labeling Buffer, 12.5ml 10X HBSS, 25ml

For Research Use Only, Not for Human

GLOBAL HEADQUARTERS

Enzo Life Sciences, Inc.
10 Executive Blvd
Farmingdale, NY 11735
USA
T 1-800-942-0430
T 1-631-694-7070
F 1-610-941-9252
E info-usa@enzolifesciences.com
www.enzolifesciences.com

EUROPE/ASIA

Enzo Life Sciences (ELS) AG
Industriestrasse 17, Postfach
CH-4415 Lausen
Switzerland
T +41/061 926 89 89
F +41/061 926 89 79
E info-ch@enzolifesciences.com
www.enzolifesciences.com




Figure 1. Fluorescence excitation (359, 460) and emission (527) spectra for the CYTO-ID[®] Green tracer dye. All spectra were determined for cell-bound dye.




Figure 2: Composite bright-field (A) and fluorescence microscopy images (B) demonstrating staining of Jurkat cells with CYTO-ID[®] Green Tracer dye. Standard FITC (Green) filter set was used to image the membrane-bound signal.

GLOBAL HEADQUARTERS

Enzo Life Sciences, Inc.
10 Executive Blvd
Farmingdale, NY 11735
USA
T 1-800-942-0430
T 1-631-694-7070
F 1-610-941-9252
E info-usa@enzolifesciences.com
www.enzolifesciences.com

EUROPE/ASIA

Enzo Life Sciences (ELS) AG
Industriestrasse 17, Postfach
CH-4415 Lausen
Switzerland
T +41/061 926 89 89
F +41/061 926 89 79
E info-ch@enzolifesciences.com
www.enzolifesciences.com




Figure 3: Flow cytometry analysis of fluorescence of mixed population of Jurkat cells over time. Jurkat cells stained with CYTO-ID[®] Green Tracer dye were mixed with an unstained population of Jurkat cells and incubated over a 120-hour period.

Product Literature References

The therapeutic response of CDDO-Me in the esophageal squamous cell carcinoma (ESCC) cells is mediated by CaMKIIa Y. Wang, et al. Am. J. Transl. Res. **8** 1695 (2016)

Is the autophagy a friend or foe in the silver nanoparticles associated radiotherapy for glioma? H. Wu, et al. Biomaterials **62** 47 (2015)

Revised 23-Dec-16

GLOBAL HEADQUARTERS

Enzo Life Sciences, Inc.
10 Executive Blvd
Farmingdale, NY 11735
USA
T 1-800-942-0430
T 1-631-694-7070
F 1-610-941-9252
E info-usa@enzolifesciences.com
www.enzolifesciences.com

EUROPE/ASIA

Enzo Life Sciences (ELS) AG
Industriestrasse 17, Postfach
CH-4415 Lausen
Switzerland
T +41/061 926 89 89
T +41/061 926 89 79
E info-ch@enzolifesciences.com
www.enzolifesciences.com