

**COLLAGEN, TYPE IV, HUMAN,
LYOPHILIZED, STERILIZED**
Catalog Number **5022**

Product Description

Type IV collagen is the primary collagen found in the extracellular basement membranes separating a variety of epithelial and endothelial cells. It is a major component of the dermal-epidermal junction where it is mostly found in the lamina densa. It is a heterotrimeric molecule containing two α 1-like and one α 2-like chains.

This Type IV collagen is isolated from human placenta and is purified using a multi-step process. The product is supplied as a sterile, lyophilized powder containing 5 mg of Type IV collagen per vial.

Type IV human collagen is typically used as a thin coating on tissue culture surfaces. This product is generally used *in vitro* as a substrate scaffold to enhance cell attachment, adherence and proliferation. Type IV collagen may be used to culture epithelial, endothelial, muscle, nerve and many other cell types. Additionally, this product is suitable for use as a substrate for collagenase assays and positive controls.

Characterization

Identity/Purity: The identity and purity of Type IV Human Collagen is qualitatively evaluated using electrophoresis (SDS-PAGE) which shows the typical banding pattern for Type IV collagen.

Storage: This product is stored at -20°C prior to solubilization and is shipped on frozen gel packs. The product is recommended to be stored at 2 to 10°C after reconstitution. Avoid repeated freeze and thaw.

Stability: The product shelf life is 24 months when stored at -20°C . The product shelf life after reconstitution is 3 months when stored at $2-10^{\circ}\text{C}$.

Cell Adherence Assay: To demonstrate the bioactivity, human dermal fibroblasts were seeded onto surfaces coated with Type IV Collagen in serum free conditions. All surfaces were blocked with a solution containing 1% BSA. Cells were then allowed to attach for one (1) hour at 37°C . The results indicate significant cell attachment bioactivity of Type IV Collagen. The control surfaces showed minimal and poor cell adherence.

Precautions and Disclaimer

This product is for R&D use only and is not intended for human or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

The raw material source for this product is human placenta. The raw material has been tested for the presence of infectious viruses (HIV 1 & 2, HBV, HCV) and found nonreactive. However, no known test method can offer complete assurance of safety. Appropriate safety and personal protective practices should be followed when handling this product.

Preparation Procedure

1. Reconstitute the 5 mg vial with 5 ml of sterile 0.25% acetic acid and mix vigorously until the contents are completely solubilized.
2. Dilute the product to desired concentration with 0.25% acetic acid. A typical final coating concentration may be 10 to $100\ \mu\text{g}/\text{cm}^2$. Testing will likely be required to determine optimal concentrations required for different cell culture systems.
3. Add appropriate amount of diluted product to culture surface.
4. Incubate at room temperature or 37°C , covered, for 1-2 hours.
5. After incubation, aspirate any excess remaining material.
6. Rinse coated surface carefully with a sterile balanced salt solution. Avoid scratching surfaces.
7. Aspirate remaining material from coated surface.
8. Coated culture vessels are now ready to use. The coated culture vessels may be stored at 2 to 10°C .