



# Product Manual

## **MULTIVIEW<sup>®</sup> (mouse-HRP/mouse-AP) IHC kit**

Catalog #: ADI-950-101



# Product Manual

## **USE FOR RESEARCH PURPOSES ONLY**

Unless otherwise specified expressly on the packaging, all products sold hereunder are intended for and may be used for research purposes only and may not be used for food, drug, cosmetic or household use or for the diagnosis or treatment of human beings. Purchase does not include any right or license to use, develop or otherwise exploit these products commercially. Any commercial use, development or exploitation of these products or development using these products without the express written authorization of Enzo Life Sciences, Inc. is strictly prohibited. Buyer assumes all risk and liability for the use and/or results obtained by the use of the products covered by this invoice whether used singularly or in combination with other products.

## **LIMITED WARRANTY; DISCLAIMER OF WARRANTIES**

These products are offered under a limited warranty. The products are guaranteed to meet all appropriate specifications described in the package insert at the time of shipment. Enzo Life Sciences' sole obligation is to replace the product to the extent of the purchasing price. All claims must be made to Enzo Life Sciences, Inc., within five (5) days of receipt of order. THIS WARRANTY IS EXPRESSLY IN LIEU OF ANY OTHER WARRANTIES OR LIABILITIES, EXPRESS OR IMPLIED, INCLUDING WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, AND NON- INFRINGEMENT OF THE PATENT OR OTHER INTELLECTUAL PROPERTY RIGHTS OF OTHERS, AND ALL SUCH WARRANTIES (AND ANY OTHER WARRANTIES IMPLIED BY LAW) ARE EXPRESSLY DISCLAIMED.

## **TRADEMARKS AND PATENTS**

Several Enzo Life Sciences products and product applications are covered by US and foreign patents and patents pending. Enzo is a trademark of Enzo Life Sciences, Inc.

**FOR RESEARCH USE ONLY.  
NOT FOR USE IN DIAGNOSTIC PROCEDURES.**

## TABLE OF CONTENTS



Please read entire booklet before proceeding with the assay.



Carefully note the handling and storage conditions of each kit component.



Please contact Enzo Life Sciences Technical Support if necessary.

Intended Use.....	2
Kit Contents .....	2
Storage .....	2
Stability .....	2
Composition.....	3
Other Materials Needed .....	3
Precautions .....	3
Recommended Staining Protocol.....	4
Staining Procedure.....	6
Contact Information.....	10



## INTENDED USE

The MULTIVIEW<sup>®</sup> (mouse-HRP/mouse-AP) IHC kit is a non-biotin, one-step detection system suitable for demonstrating antigens in formalin-fixed paraffin-embedded tissues and cryostat sections. It may also be used with blood smears, cytosmears, and cell preparations.

The MULTIVIEW<sup>®</sup> (mouse-HRP/mouse-AP) IHC kit has been developed by directly labeling anti-mouse immunoglobulins with HRP and AP enzymes using a proprietary tandem hyperlabelling technology. This ensures consistent and reproducible immunodetection of 2 mouse monoclonal antibodies. Nuclear, cytoplasmic, and membrane antigens in different types of tissues can be readily detected.

The MULTIVIEW<sup>®</sup> (mouse-HRP/mouse-AP) IHC kit allows the user to detect two different mouse antibodies on the same section via sequential staining and chromogen development. The reagents can be used for manual staining or with automated staining instruments and are specifically designed for multiplex immunohistochemical staining assays.

## KIT CONTENTS

- |  |                                |
|--|--------------------------------|
| 1. Peroxidase Block  | 10 mL                          |
| 2. POLYVIEW <sup>®</sup> IHC reagent (mouse-HRP) (Prod. no. ADI-950-112)                                   | 10 mL                          |
| 3. POLYVIEW <sup>®</sup> IHC reagent (mouse-AP) (Prod. no. ADI-950-110)                                    | 10 mL                          |
| 4. HIGHDEF <sup>®</sup> IHC chromogen substrate (DAB, stable) (Prod. no. ADI-950-212) buffer and chromogen | 15 mL / 1 mL<br>15 mL / 0.5 mL |
| 5. HIGHDEF <sup>®</sup> red IHC chromogen (AP) (Prod. no. ADI-950-140) buffer & chromogen                  |                                |
| 6. Empty Mixing Bottles  | 2 bottles                      |

## STORAGE

Store kit and reagents at 2-8°C in the dark. Do not freeze.

## STABILITY

12-24 months (see expiration date on reagent bottles).



# Product Manual

## COMPOSITION

Peroxidase Block, The POLYVIEW<sup>®</sup> detection reagents, Mouse HRP and Mouse AP, are provided in ready to use format. HIGHDEF<sup>®</sup> IHC chromogen substrate (DAB, stable), and HIGHDEF<sup>®</sup> red IHC chromogen (AP) solutions require preparation before use.

Refer to "Preparation" sections for this information.

## OTHER MATERIALS NEEDED

- Xylene or dewaxing reagents
- Absolute ethanol
- Distilled or deionized water
- IHC wash buffer (Prod. no. ADI-950-235)
- Pre-Blocking solution: IHC background blocker (Prod. no. ADI-950-230, optional)
- IHC diluent (primary Ab) (Prod. no. ADI-950-244)
- Counterstain
- Mounting Medium: HIGHDEF<sup>®</sup> mount (Prod. no. ADI-950-261)

## PRECAUTIONS

1. Wear appropriate personal protective apparel. Avoid contact with clothes and exposed skin. In case of accidental skin exposure, flush with water immediately. Consult a physician if required.
2. Interpretation of the results is the sole responsibility of the user.

## RECOMMENDED STAINING PROTOCOL

1. Paraffin-embedded tissue sections must be deparaffinized with xylene or dewaxing agent and rehydrated with a graded series of ethanol and water washes before staining. Follow the standard dewaxing and rehydration protocol used in your lab.
2. The investigator needs to optimize the dilution and incubation times for primary antibodies.
3. Each immunostaining run should include known positive and negative controls to assure proper functioning of the staining system and aid in valid interpretation of the results.

### Typical controls:

*Positive Control:* A tissue known to contain the desired antigen, which has yielded positive staining in the past.

*Negative Controls:*

### Reagent Controls

- A. Substitute normal non-immune serum from the same host animal as the primary antibody (e.g. if using mouse monoclonal primary antibodies, use mouse non-immune serum).
- B. Substitute matching host species isotype control for primary antibody
- C. Use antigen-adsorbed primary antibody (i.e. antibody reagent which has been adsorbed with the target antigen to remove specific antibody)

Tissue control – A tissue known to not contain the desired antigen.

4. Consult the primary antibody supplier for recommended antigen recovery treatments. Perform epitope recovery pretreatments before starting staining procedures.
5. Once the slide treatment has been started, DO NOT let tissues or specimens dry. This can cause undesirable background or artifacts.



## **Preparation of HIGHDEF® IHC chromogen substrate (DAB, stable) Working Solution**

1. Transfer 1mL of HIGHDEF® IHC chromogen substrate (DAB, stable) buffer to a tube or mixing bottle.
2. Add 1 drop (approximately 20  $\mu$ L) of HIGHDEF® IHC chromogen substrate (DAB, stable) chromogen to the buffer and mix thoroughly.
3. HIGHDEF® IHC chromogen substrate (DAB, stable) working solution is stable for 2 weeks refrigerated at 2-8°C.
4. Working solution volume can be scaled up using the same ratio of buffer to chromogen.
5. Dispose of unused HIGHDEF® IHC chromogen substrate (DAB, stable) working solution in appropriate waste stream according to local, state or federal regulations.

## **Preparation of HIGHDEF® red IHC chromogen (AP) Substrate Working Solution**

1. Transfer 1 mL of HIGHDEF® red IHC chromogen (AP) Buffer to a tube or mixing bottle.
2. Add 1 drop (approximately 20  $\mu$ L) of HIGHDEF® red IHC chromogen (AP) chromogen to buffer and mix thoroughly.
3. HIGHDEF® red IHC chromogen (AP) substrate working solution should be used within 20-30 minutes of preparation. Discard any solution not used within this period.
4. Working solution volume can be scaled up using the same ratio of buffer to chromogen.
5. Dispose of unused HIGHDEF® red IHC chromogen (AP) substrate working solution in appropriate waste stream according to local, state, and federal regulations.

## STAINING PROCEDURE

		INCUBATION TIME
<b>1. Perform antigen retrieval (if necessary)</b>	Recommended solutions: IHC enzyme antigen retrieval reagent (Prod. No. ADI-950-280), Antigen Retrieval Reagent, pH 6 (10X) (Prod. No. ENZ-ACC112), or Antigen Retrieval Reagent, pH 9 (10X) (Prod. No. ENZ-ACC113).	<b>5-25 min.</b>
<b>2. Block endogenous HRP activity</b>	Incubate with Peroxidase Block Wash slides with 3 changes of IHC wash buffer.	<b>5– 10 min.</b> <b>3 x 1 min.</b>
<b>3. Apply first primary mouse antibody</b>	Apply mouse primary antibody that will be used with HRP. Wash slides with 3 changes of IHC wash buffer.	<b>30 min.</b> <b>3 x 1 min.</b>
<b>4. Apply Mouse HRP Solution</b>	Apply POLYVIEW® IHC reagent (Mouse HRP) solution. Wash slides with 3 changes of IHC wash buffer.	<b>20 min.</b> <b>3 x 1 min.</b>
<b>5. Develop with HIGHDEF® substrate (DAB, stable)</b>	Apply the HIGHDEF® IHC chromogen substrate (DAB, stable) solution and incubate. Wash slides with 3 changes of IHC wash buffer.	<b>5 min.</b> <b>3 x 1 min.</b>



<b>6. Apply second primary mouse antibody</b>	<p>Apply the second mouse primary antibody that will be used with AP.</p> <p>Wash slides with 3 changes of IHC wash buffer.</p>	<p><b>30 min.</b></p> <p><b>3 x 1 min.</b></p>
<b>7. Apply Mouse AP Solution</b>	<p>Apply POLYVIEW<sup>®</sup> IHC reagent (Mouse-AP) solution.</p> <p><b>Wash slides with 3 changes of IHC wash buffer.</b></p>	<p><b>20 min.</b></p> <p><b>3 x 1 min.</b></p>
<b>8. Develop with HIGHDEF<sup>®</sup> Red (AP)</b>	<p>Apply the HIGHDEF<sup>®</sup> red IHC chromogen (AP) solution and incubate.</p> <p><b>Wash slides with 3 changes of IHC wash buffer.</b></p>	<p><b>15 min.</b></p> <p><b>3 x 1 min.</b></p>

<b>9. Counterstain</b>	<p>Incubate tissue with Counterstain (e.g. Hematoxylin), according to manufacturer's recommendation or standard laboratory protocol.</p> <p><b>Wash slides with water 3 times, followed by 1 time in IHC wash buffer, then 1 time in water.</b></p>	<p><b>1-5 min.</b></p> <p><b>5 x 1 min.</b></p>
<b>10. Dehydrate &amp; Coverslip</b>	<p>Dehydrate tissues through graded ethanol series, followed by xylene series.</p> <p>Apply coverslips with permanent mounting medium.</p>	

**NOTE:** The use of an antibody removal buffer between the first layer and the second is optional. Although, sufficient washing, as indicated in the staining protocol (3 x 1 minute washing with IHC wash buffer), will remove excess antibody, detection reagent, or chromogen making the staining of each antigen of interest specific.



# Product Manual

**NOTES**



# Product Manual

## **GLOBAL HEADQUARTERS**

Enzo Life Sciences Inc.  
10 Executive Boulevard  
Farmingdale, NY 11735  
Toll-Free: 1.800.942.0430  
Phone: 631.694.7070  
Fax: 631.694.7501  
info-usa@enzolifesciences.com

## **EUROPE/ASIA**

Enzo Life Sciences (ELS) AG  
Industriestrasse 17  
CH-4415 Lausen  
Switzerland  
Phone: +41/0 61 926 89 89  
Fax: +41/0 61 926 89 79  
info-ch@enzolifesciences.com

For local distributors and detailed product information visit us online:

[www.enzolifesciences.com](http://www.enzolifesciences.com)