



> **ImmunoSet® HSP25 (rodent),
ELISA development set**
Catalog # ADI-960-075
Reagents for 5 x 96-Well EIA Kits

This ImmunoSet contains the basic components for the development of an Hsp25 immunometric enzyme immunoassay (EIA). Each kit contains sufficient reagents for five 96-well plates. This kit has been validated for use with cell lysates. Additional sample types will require validation by the user. Visit www.enzolifesciences.com for tips and frequently asked questions.

Introduction

Rodent Hsp25, human Hsp27, and α B-crystallin are part of a diverse family of small heat shock proteins (sHsps) that range in size from 14-45 kDa¹. sHsps are constitutively expressed, and upregulated in response to stress such as heat shock. They function as chaperone-like proteins by binding unfolded polypeptides and preventing uncontrolled protein aggregation. Data indicates that Hsp25 is a dynamic tetramer of tetramers with a unique ability to refold and reassemble into its active quaternary structure after denaturation². Hsp25 helps facilitate the glutathione-redox cycle by enhancing glutathione utilization and maintaining the cellular glutathione pool in favor of the reduced states³.

References:

1. Narberhaus, F. (2002) *Microbiol Mol Biol Rev.* **66**, 64-93.
2. Ehrnsperger, M., *et al.* (1999) *J Biol Chem.* **274**, 14867-14874.
3. Baek, S.H., *et al.* (2000) *J Cell Physiol* **183**, 100-107.

Materials Provided

1. Hsp25 Capture Antibody
One vial containing 312.5 μ g lyophilized Hsp25 monoclonal antibody, Cat. #80-1942
2. Hsp25 Standard
One vial containing 312.5 ng lyophilized recombinant Hsp25 protein, Cat. #80-1943
3. Hsp25 Detection Antibody
One vial containing 31.25 μ g lyophilized Hsp25 biotinylated polyclonal antibody, Cat. #80-1944
4. SA-HRP
One vial containing 12.5 μ g lyophilized streptavidin conjugated to horseradish peroxidase, Cat. #80-1896

Materials Needed but not Supplied

1. RIPA Cell Lysis Buffer, Cat. #80-1284, or similar
2. 96-well high-binding polystyrene microtiter plates, Cat. #80-1930, or similar
3. Precision pipets
4. Microplate reader capable of reading at 450 nm
5. Microplate shaker
6. Phosphate buffered saline (PBS)[†]
7. Tween[®]-20*[†]
8. Bovine Serum Albumin (BSA)[†]
9. 3,3',5,5' tetramethylbenzidine (TMB) solution, Cat. #80-1805 or similar[†]
10. 1N hydrochloric acid, such as Stop Solution 2, Cat. #80-1804[†]
11. Sucrose
[†]ImmunoSet Buffer Pack, Cat. #ADI-950-003
^{*}Tween is a registered trademark of ICL Americas

Buffer Formulations

1. Coating Buffer
10 mM sodium phosphate, 15 mM NaCl, pH 7.4
2. Blocking Buffer
10 mM sodium phosphate, 15 mM NaCl, 1.0% BSA, 1% sucrose, pH 7.4
3. Assay Buffer
100 mM sodium phosphate, 150 mM NaCl, 1.0% BSA, 0.1% Tween-20, pH 7.4
4. Wash Buffer
10 mM sodium phosphate, 15 mM NaCl, 0.1% Tween-20, pH 7.4

Plate Coating

1. Reconstitute Hsp25 Capture Antibody with 250 μ L deionized water for a 250x stock. Use immediately, or make aliquots and freeze at -20°C for up to 3 months. For prolonged storage, aliquot and freeze at -80°C. Avoid repeated freeze/thaw cycles.
2. Dilute the stock 1:250 in Coating Buffer. Immediately dispense into 96-well microtiter plates using 100 μ L of the diluted capture antibody per well. Seal the plate and incubate overnight at room temperature.
3. Aspirate each well to remove coating solution. Immediately add 200 μ L Blocking Buffer per well. Seal the plate and incubate for at least 1 hour.
4. Aspirate each well to remove blocking solution. Plates may be used immediately or dried and stored with desiccant at 4°C.

Reagent Preparation

1. Recombinant Hsp25 Standard
Reconstitute vial contents with 250 μ L deionized water for a 1250 ng/mL (50x) stock. Aliquot and store at -20°C for up to 3 months. For prolonged storage, aliquot and freeze at -80°C. Avoid repeated freeze/thaw cycles.
The recommended standard curve range is 25 ng/mL to 0.78 ng/mL, using 2-fold serial dilutions in Assay Buffer. Do not store diluted standard.
2. Hsp25 Detection Antibody
Reconstitute vial contents with 250 μ L deionized water for a 250x stock. Aliquot and store at -20°C for up to 3 months. For prolonged storage, aliquot and freeze at -80°C. Avoid repeated freeze/thaw cycles.
Dilute the stock 1:250 in Assay Buffer for a working solution. Do not store diluted antibody.
3. SA-HRP
Reconstitute vial contents with 250 μ L deionized water for a 1500x stock. Store at 4°C for up to 3 months. For prolonged storage, aliquot and freeze at -20°C. Avoid repeated freeze/thaw cycles.
Dilute the stock 1:1500 in Assay Buffer for a working solution. Do not store diluted conjugate.

Assay Procedure

1. Pipet 100 μ L of Assay Buffer into the control (0 ng/mL standard) wells.
2. Pipet 100 μ L of standards and samples, prepared in Assay Buffer, to the bottom of the appropriate wells.
3. Seal the plate. Incubate for 1 hour on a plate shaker at room temperature.
4. Empty the contents of the wells and wash by adding 400 μ L of Wash Buffer to every well. Repeat 3 more times for a total of 4 washes. After the final wash, empty or aspirate the wells and firmly tap the plate on a lint free paper towel to remove any remaining wash buffer.
5. Pipet 100 μ L of the diluted detection antibody into each well, except the blank.
6. Seal the plate. Incubate for 1 hour on a plate shaker at room temperature.
7. Wash as above (Step 4).
8. Add 100 μ L of the diluted conjugate to each well except the blank.
9. Seal the plate. Incubate for 30 minutes on a plate shaker at room temperature.
10. Wash as above (Step 4).
11. Pipet 100 μ L of TMB solution into each well.
12. Seal the plate. Incubate for 30 minutes on a plate shaker at room temperature.
13. Pipet 100 μ L 1N HCl into each well.
14. After blanking the plate reader against the substrate, read optical density at 450 nm. If the plate reader is not capable of adjusting for the blank, manually subtract the mean OD of the substrate blank from all readings.

