



> ImmunoSet® Osteopontin (human), ELISA development set
Catalog # ADI-960-142
Reagents for 5 x 96-Well EIA Kits

This ImmunoSet contains the basic components for the development of an OPN (human) immunometric enzyme immunoassay (EIA). Each kit contains sufficient reagents for ten 96-well plates.

This kit has been validated for use with cell culture supernatants, EDTA and citrate plasma, and milk. Additional sample types will require validation by the user.

Visit www.enzolifesciences.com for tips and frequently asked questions.

Introduction

Osteopontin (OPN) is an acidic extracellular matrix cell adhesion protein abundant in bone matrix, plasma, urine, and milk, which directs a variety of processes including tissue remodeling, inflammation, and cell survival¹. Plasma OPN has been shown to be a positive indicator of breast, colon, and lung cancers as well as metastatic carcinomas^{2,3}. The presence of OPN in a variety of tumors is strongly correlated to pathological stage, suggesting its critical role in tumor invasiveness, progression and metastasis^{4,5}. OPN is found in atherosclerotic plaques and may drive a number of diabetic vascular pathologies⁶.

1. Denhardt, D.T., *et al.* (2001) *J Clin Invest.* **107**, 1055-1061.
2. Agrawal, D., *et al.* (2002) *J Natl Cancer Inst.* **94**, 513-521.
3. Hotte, S.B., *et al.* (2002) *Cancer* **95**, 506-512.

4. Coppola, D., *et al.* (2004) *Clin Cancer Res.* **10**, 184-190.
5. Tuck, A.B., *et al.* (1997) *Arch Pathol Med.* **121**, 578-584.
6. Takemoto, M., *et al.* (2000) *Arterioscler Thromb Vasc Biol.* **20**, 624-628.

Materials Provided

1. OPN (human) Capture Antibody
Two vials containing 344 µg each lyophilized OPN (human) monoclonal antibody, Cat. #80-1956
2. OPN (human) Standard
Two vials containing 12.5 ng each lyophilized recombinant OPN (human) protein, Cat. #80-1957
3. OPN (human) Detection Antibody
Two vials containing 9.4 µg each lyophilized OPN (human) polyclonal antibody, Cat. #80-1958
4. SA-HRP
One vial containing 12.5 µg lyophilized streptavidin conjugated to horseradish peroxidase, Cat. #80-1896

Materials Needed but not Supplied

1. 96-well high-binding polystyrene microtiter plates, Cat. #80-1930, or similar
2. Precision pipets
3. Microplate reader capable of reading at 450 nm
4. Phosphate buffered saline (PBS)[†]
5. Tween®-20*[†]
6. Bovine Serum Albumin (BSA)[†]
7. 3,3',5,5' tetramethylbenzidine (TMB) solution, Cat. #80-1805 or similar[†]
8. 1N hydrochloric acid, such as Stop Solution 2, Cat. #80-1804[†]
9. Sucrose

[†]ImmunoSet Buffer Pack, Cat. #ADI-950-003

*Tween is a registered trademark of ICL Americas

Buffer Formulations

1. Coating Buffer
10 mM sodium phosphate, 15 mM NaCl, pH 7.4

2. Blocking Buffer
10 mM sodium phosphate, 15 mM NaCl, 1.0% BSA, 1.0% Sucrose, pH 7.4
3. Assay Buffer
10 mM sodium phosphate, 15 mM NaCl, 1.0% BSA, 0.1% Tween-20, pH 7.4
4. Wash Buffer
10 mM sodium phosphate, 15 mM NaCl, 0.1% Tween-20, pH 7.4

Plate Coating

1. Reconstitute OPN (human) Capture Antibody with 250 µL deionized water for a 250x stock. Use immediately, or make aliquots and freeze at -20°C for up to 3 months. For prolonged storage, aliquot and freeze at -80°C. Avoid repeated freeze/thaw cycles.
2. Dilute the stock 1:250 in Coating Buffer. Immediately dispense into 96-well microtiter plates using 100 µL of the diluted capture antibody per well. Seal the plate and incubate overnight at room temperature.
3. Aspirate each well to remove coating solution. Immediately add 200 µL Blocking Buffer per well. Seal the plate and incubate for at least 1 hour.
4. Aspirate each well to remove blocking solution. Plates may be used immediately or dried and stored with desiccant at 4°C.

Reagent Preparation

1. Recombinant OPN (human) Standard
Reconstitute vial contents with 250 µL deionized water for a 50 ng/mL (50x) stock. Aliquot and store at -20°C for up to 3 months. For prolonged storage, aliquot and freeze at -80°C. Avoid repeated freeze/thaw cycles.
The recommended standard curve range is 1.0 ng/mL to 0.03 ng/mL, using 2-fold serial dilutions in Assay Buffer. Do not store diluted standard.
2. OPN (human) Detection Antibody
Reconstitute vial contents with 250 µL deionized water for a 250x stock. Store at 4°C for up to 3 months. For prolonged storage, aliquot and freeze at -20°C. Avoid repeated freeze/thaw cycles.

Dilute the stock 1:250 in Assay Buffer for a working solution. Do not store diluted antibody.

3. SA-HRP
Reconstitute vial contents with 250 µL deionized water for a 667x stock. Store at 4°C for up to 3 months. For prolonged storage, aliquot and freeze at -20°C. Avoid repeated freeze/thaw cycles.
Dilute the stock 1:667 in Assay Buffer for a working solution. Do not store diluted conjugate.


Assay Procedure

1. Pipet 100 µL of Assay Buffer into the control (0 ng/mL standard) wells.
2. Pipet 100 µL of standards and samples, prepared in Assay Buffer, to the bottom of the appropriate wells.
3. Seal the plate. Incubate on a plate shaker for 1 hour at room temperature.
4. Empty the contents of the wells and wash by adding 400 µL of Wash Buffer to every well. Repeat 3 more times for a total of 4 washes. After the final wash, empty or aspirate the wells and firmly tap the plate on a lint free paper towel to remove any remaining wash buffer.
5. Pipet 100 µL of the diluted detection antibody into each well, except the blank.
6. Seal the plate. Incubate on a plate shaker for 1 hour at room temperature.
7. Wash as above (Step 4).
8. Add 100 µL of the diluted conjugate to each well except the blank.
9. Seal the plate. Incubate on a plate shaker for 30 minutes at room temperature.
10. Wash as above (Step 4).
11. Pipet 100 µL of TMB solution into each well.
12. Seal the plate. Incubate on a plate shaker for 30 minutes at room temperature.
13. Pipet 100 µL 1N HCl into each well.
14. After blanking the plate reader against the substrate, read optical density at 450 nm. If the plate reader is not capable of adjusting for the blank, manually subtract the mean OD of the substrate blank from all readings.

Assay Performance

Typical Data

The results shown below are for illustration only and should not be used to interpret results from another assay.



Sensitivity

The sensitivity, or limit of detection, of this assay is 0.016 ng/mL. It was determined by interpolation at 2 standard deviations above the mean signal at background, using data from 7 standard curves.

Specificity

This assay detects OPN in cell culture supernatants, Na-Citrate and EDTA plasmas, and milk of human origin. Cross reactivity with mouse and bovine OPN is less than 0.02%.

Dilutional Linearity

To determine possible interference from the sample matrix, the indicated sample types were serially diluted into assay buffer. The concentrations of OPN were measured in the assay, and the results were analyzed to determine the range over which a linear response was obtained. From these data, the minimum recommended dilution (MRD) is as follows: culture supernatants (1:32), plasma (1:64), and milk (1:8,000).


Dilution Factor*	MDA-MB-435 CS	Citrate Plasma	Milk
1:8	---	---	100%
1:16	---	---	130%
1:32	100%	---	107%
1:64	93%	100%	91%
1:128	81%	133%	---
1:256	83%	137%	---
1:512	---	103%	---

CS: Culture Supernatant

* x1000 for milk.

Parallelism

Dose-response curves from cell lysates diluted into assay buffer (using the MRD) were compared to the recombinant OPN standard curve. A parallel response indicates the recombinant standard effectively mimics the native protein.



Calculation of Results

Several options are available for the calculation of the relative levels of OPN in samples. We recommend that the data be handled by an immunoassay software package utilizing a 4 parameter logistic curve-fitting program. For accuracy, please ensure that sample values fall within the standard range.

Accessory Reagent List		
Reagent	Quantity	Cat. #
ImmunoSet® Buffer Pack	1 each of the following products: 80-1927, 80-1928, 80-1929, 80-1805, 80-1804	ADI-950-003
ImmunoSet® Plate Pack	5 96-well clear microtiter plates & 5 plate sealers	80-1930
PBS Concentrate	120 mL	80-1927
BSA Solution (10%)	50 mL	80-1928
Tween-20 Solution (10%)	30 mL	80-1929
Wash Buffer Concentrate	100 mL	80-1287
SA-HRP	12.5 µg/vial	80-1896

Storage

Store all components at 4°C. See page 3 for storage of reconstituted material.

Tips & Troubleshooting

- ✓ If buffers other than those recommended are used in the assay, the end-user must determine the appropriate dilution and assay validation.
- ✓ Pipet the reagents to the sides of the wells to avoid possible contamination.
- ✓ Pre-rinse each pipet tip with reagent. Use fresh pipet tips for each sample, standard, and reagent.
- ✓ Insufficient washing or residual wash buffer in the wells may cause variation in assay results.
- ✓ Bring all reagents to room temperature for at least 30 minutes prior to opening.
- ✓ All standards, controls, and samples should be assayed in duplicate.

Limited Warranty

Enzo Life Sciences International, Inc. makes no warranty of any kind, expressed or implied, which extends beyond the description of the product in this brochure, except that the material will meet our specifications at the time of delivery. Enzo Life Sciences International, Inc. makes no guarantee of results and assumes no liability for injuries, damages or penalties resulting from product use, since the conditions of handling and use are beyond our control.