

**PERKINELMER™ LIFE SCIENCES, INC.**



**OLIGONUCLEOTIDE 5' END  
LABELING SYSTEM**

**Catalog Number NEP101**

**For Laboratory Use Only  
Caution: Research chemicals for research purposes only.**



Contains reagents for 5' end labeling oligonucleotides.

**Note: Tracer must be ordered separately.**

## I. INTRODUCTION

The techniques for end labeling oligonucleotides with radioisotopes have allowed for rapid developments in nucleic acid probe technology. Oligonucleotide probes can be custom made based on sequence information of the target DNA or RNA in several hours on a DNA synthesizer. This eliminates the usual cumbersome and time consuming steps that are involved in the cloning and isolation of restriction fragments to be used as hybridization probes. Another advantage of oligonucleotide probes is that they can be designed to detect single base changes (mutations) in a gene<sup>1-5</sup>. This is an especially useful property, by adjusting the conditions of the hybridization and washes (high stringency), only probes with a perfect match can be made to hybridize. Several techniques have been developed which allow single base changes to be detected using these DNA or RNA probes<sup>6,7</sup>.

Synthetic oligonucleotides are prepared with free 5' and 3' hydroxyl groups. This means there is no need for pretreatment with bacterial alkaline or calf intestinal phosphatase. End labeling of oligonucleotides is also much simpler since there is no need to modify reaction conditions depending upon whether the fragment has 5' or 3' overhangs or blunt ends as with restriction fragments<sup>8</sup>. Oligonucleotides can be labeled at either the 3' or 5' end. Using polynucleotide kinase and ATP <sup>32</sup>P (BLU/NEG502Z) the 5' end is labeled. Using terminal transferase and deoxynucleotide 5' triphosphates the 3' end is labeled. The isotope of choice has been <sup>32</sup>P, however <sup>35</sup>S and <sup>125</sup>I have been used successfully<sup>9</sup>. Oligonucleotides labeled with <sup>32</sup>P are reported to be stable for 1-2 weeks and those labeled with <sup>35</sup>S last approximately 6 times longer. The stability of <sup>125</sup>I probes is currently under evaluation. The <sup>35</sup>S or <sup>125</sup>I labeled oligonucleotides are especially useful when high resolution (as in *in situ* hybridization) or long probe stability is needed.

## II. PRINCIPLE OF OLIGONUCLEOTIDE 5' END LABELING

Polynucleotide kinase, simultaneously identified by Richardson<sup>10</sup>, and Novagrodsky and Hurwitz<sup>11</sup> in T-phage infected *E. coli*, catalyzes the transfer of the gamma phosphate of ATP to the 5' hydroxyl terminus of DNA or oligonucleotide molecules. Oligonucleotides can be labeled with radioactivity if the gamma phosphate of ATP contains phosphorous-32. Labeling the 5' ends with <sup>35</sup>S utilizing ATP (gamma thiophosphate), [<sup>35</sup>S] is not recommended due to the slow rate of the enzyme with ATP $\alpha$ S as substrate<sup>12</sup>. Complete phosphorylation of the 5' hydroxyls is achieved using low concentrations of ATP, usually a two to five-fold excess over the concentration of 5' hydroxyl ends.

### III. EXPLANATION OF THE SYSTEM

PerkinElmer Life Sciences, Inc. End Labeling System is an adaptation of the procedure described by Maniatis<sup>13</sup>, which is optimized to obtain maximal labeling of the oligonucleotide. The system includes buffers, enzyme, a control oligonucleotide and the other reagents for labeling the 5' termini of oligonucleotides. To check that the system is functioning properly, a 17 mer oligonucleotide has been included as a control. The T<sub>4</sub>-polynucleotide kinase has been purified to homogeneity, as demonstrated by SDS. gel electrophoresis, and has been shown to be free of any detectable contaminating ribonuclease, deoxyribonuclease, or phosphatase activities.\*

T<sub>4</sub>-Polynucleotide Kinase, NEE101.

### IV. THE FOLLOWING REAGENTS ARE SUPPLIED IN THE KIT

Store kit components at temperature indicated in this section, in general, all enzymes and reagents at -20°C. Enough of each reagent is supplied with this kit for minimally 10 labeling reactions.

#### T<sub>4</sub> Polynucleotide Kinase (NEE101)

One vial containing 30 ul (30 units\*) of polynucleotide kinase (EC2.7.1.78) from T<sub>4</sub>-infected *E. coli* is provided already diluted to the proper concentration for use (1000 units/ml). The enzyme has been purified to homogeneity and is stocked in Tris-HCl buffer, pH7.6, containing ATP, DTT, KCl, EDTA, and 50% glycerol. The material is stable for at least one month when stored at -20°C.

#### Phosphorylation Buffer

One vial containing 0.15 ml of a 10x buffer is supplied for the phosphorylation reaction. The buffer is a solution of Tris-HCl, pH7.6, MgCl<sub>2</sub>, DTT, and spermidine. The material is stable for at least one month when stored at -20°C.

#### Control Oligonucleotide

Approximately 50 pmole of a 17 mer control oligonucleotide is supplied in 50 µl of deionized water. The preparation is supplied ready to use as a control to check the system according to the protocol provided. The material is stable at least one month when stored at -20°C.

#### Deionized Water

One vial containing 2 mL deionized water is supplied. This component is provided to minimize interfering substances (e.g. metals, ions, and

enzymes) present in some water supplies, which might inhibit the enzymatic reactions.

## **V. SUPPLIES AND EQUIPMENT NEEDED**

Vortex mixer  
Pipetting equipment  
Water bath for 37°C incubations  
Plastic gloves  
Polypropylene microcentrifuge tubes, or siliconized glass  
Brinkman PEI cellulose thin layer strips (optional)  
3 or 5 cc disposable syringes  
Microcentrifuge  
Methanol  
20% ethanol  
Water  
Isotope

## **VI. REACTION PROTOCOLS FOR LABELING 5' END LABELING**

Described below is the recommended procedure for 5' end labeling oligonucleotides with phosphorous-32. The 17mer control oligonucleotide supplied with the kit is used as a control for the labeling procedure.

### **A. Preparation of Oligonucleotide for End Labeling**

Prior to the end labeling, it is recommended that the oligonucleotide be purified to remove any contaminating salts, reagents, organic solvents, or proteins which might affect the reaction. The oligonucleotide may be purified by a wide variety of standard procedures, including HPLC or purify by method of choice. After purification, the oligonucleotide should be taken to dryness (nitrogen stream or under vacuum).

### **B. 5' End Labeling Reaction**

1. The procedure is written for 10 pmoles of oligonucleotide. The material may be in any volume up to 25  $\mu$ l of water or buffer, pH 7.6.
2. Add 5  $\mu$ l of phosphorylation buffer (10x) to the tube containing 10 pmoles of oligonucleotide.
3. Add deionized water so that the final reaction volume will be 50  $\mu$ l.
4. Add tracer to the reaction tube. The number of

picomoles of ATP, [ $\gamma$ - $^{32}\text{P}$ ], should be at least 2.0 times the number of picomoles of 5' ends. For labeling 10 pmoles of oligonucleotides, 20 pmoles of  $^{32}\text{P}$ -ATP is needed or approximately 10  $\mu\text{l}$  of ATP, [ $\gamma$ - $^{32}\text{P}$ ], (BLU/NEG502Z), solution is recommended. Mix the contents well and centrifuge briefly.

5. Add 3  $\mu\text{l}$ , (3 units) of polynucleotide kinase solution. Mix the solution gently. Spin down briefly in a microcentrifuge.
6. Incubate at 37°C for 30 minutes.
7. The reaction may be followed analytically as follows. Spot small aliquots (0.1  $\mu\text{l}$  volume, the smaller volume the better) on a PEI-Cellulose thin layer strip (Brinkman). The TLC is developed in 1M  $\text{KH}_2\text{PO}_4$  buffer for 20 minutes. The buffer in the chamber should not reach the application spot. Labeled oligonucleotide remains at the origin and ATP, [ $^{32}\text{P}$ ] migrates with an Rf of approximately 0.6-0.7. The TLC strip is either autoradiographed or scanned with a GM-tube.
8. To stop the reaction add 5  $\mu\text{l}$  of 250 mM EDTA and cool to 4°C on ice.

#### **5' END LABELING REACTION SUMMARY**

	<u>Vol. (<math>\mu\text{L}</math>)</u>
Oligonucleotide (10 pmol)*	X = 5-25
10x Phosphorylation Buffer	5
ATP, [ $\gamma$ - $^{32}\text{P}$ ], (100 uci)	10
Water	32 - X 22-2
Mix and spin down	
Add Polynucleotide Kinase	<u>3</u>
	50 $\mu\text{l}$

Incubate at 37°C for 30 minutes. Stop the reaction with 5  $\mu\text{l}$  of 0.1M EDTA pH 8.0. Place on ice.

\*1 O.D. = 20  $\mu\text{g}/\text{mL}$ , average molecular weight dNMP = 325<sup>13</sup>.  
Purify by method of choice.

## VII. REFERENCES

1. Wallace, R.B., Shaffer, J., Murphy, R.F., Bonner, J., Hirose, T., and Itakura, K. (1979) *Nuc. Acid Res.* 6, 3543-3553. Hybridization of synthetic Oligodeoxyribonucleotides to X174 DNA: The effect of single base pair mismatch.
2. Conner, B.J., Reyes, A.A., Marin, C., Itakura, K., Teplitz, R.L., and Wallace, R.B. (1983) *PNAS* 80, 278-282. Detection of sickle cell  $\beta^S$ -globin allele by hybridization with synthetic oligonucleotides.
3. Wallace, R.B., Johnson, M.J., Hirose, T., Miyake, T., Kawashima, E. H. and Itakura, K. (1981) *Nuc. Acid Res.* 9, 879-895. The use of synthetic oligonucleotides as hybridization probes II. Hybridization of oligonucleotides of mixed sequences to rabbit  $\gamma$ -globin DNA.
4. Habden, A.N., Read, M.J., Dykes, C.W., and Harford, S. (1985) *Anal. Biochem.* 144, 75-78. M13 Clones Carrying Point Mutations: Identification by Solution Hybridization.
5. Bos, J.L., Verlaain-deVries, M., Jansen, A.M., Veeneman, G.H., van Boon, J.H., and van der Eb, A.J. (1984) *Nuc. Acid Res.* 12, 9155-9163. Three different, mutations in codon 61 of the hman N-ras gene detected by oligonucleotide hybridization.
6. Wood, W.I., Gitschier, J., Lasky, L.A., and Lawn, R.M. (1985) *Proc. Natl. Acad. Sci., USA*, 82, 1585-1588. Base composition-independent oligonucleotide screening of highly complex gene libraries.
7. Winter, E., Yamamoto, F., Almoguera, C., and Perudo, M. (1985) *PNAS* 82, 7575-7579. A method to detect and characterize point mutations in transcribed genes: Amplification and over expression of the mutant c-Ki-ras allele in human tumor cells.
8. Maxam, A.M., and Gilbert, W. (1980) *Methods in Enzymology* 65, 499-560. Sequencing End-labeled DNA with Base-Specific Chemical Cleavages.
9. Collins, M.L. and Hunsaker, W.R. (1965) *Anal. Biochem.* 151, 211-224. Improved Hybridization Assays Employing Tailed Oligonucleotide Probes: A direct comparison with 5'-End-Labeled Oligonucleotide Probes and Nick-Translated



Plasmid Probes.

10. Richardson, C.C. (1965) PNAS 54, 158-165. Phosphorylation of nucleic acid by an enzyme from T<sub>4</sub> bacteriophage-infected Echerichia coli.
11. Njovagrodsky, A. and Hurwitz, J. (1966) J. Biol. Chem. 241, 2923-2932. The enzymic phosphorylation of ribonucleic acid and deoxyribconucleic acid I. Phosphorylation at 5'-hydroxyl termini.
12. Beltz, W.R and O'Brien, K.J. (1981) Fed.40, 1849. End Labeling of DNA restriction fragments using [35S] adenosine 5'-(gamma thio) triphosphate followed by mercury affinity chromatography.
13. Maniatis, T., Fritsch, E.F., and Sambrook, J. (1982) Molecular Cloning: A Laboratory Manual. Cold Spring Habor Lab, Cold Spring Harbor, NY.
14. Johnson, M.T., Read, B.A, Monko, A.M., Pappas, G., and Johnson, B.A. (1986) Biotedmiques 4, 64-70. A Convenient, New Method for Desalting, Deproteinizing, and Concentrating DNA or RNA.

Manufactured by:

PerkinElmer Life Sciences, Inc.  
549 Albany Street  
Boston, MA USA 02118

Toll-Free 800 551-2121  
International 617 482-9595







Worldwide Headquarters: PerkinElmer Life Sciences, Inc.  
549 Albany Street, Boston, MA 02118-2512 USA (800) 551-2121  
European Headquarters: PerkinElmer Life Sciences, Imperiastraat 8, B-1930 Zaventem Belgium  
Technical Support: in Europe: techsupport.europe@perkinelmer.com in US and Rest of World:  
techsupport@perkinelmer.com Belgium: Tel: 0800 94 540 Fax: 0800 93 040  
France: Tel: 0800 90 77 62 Fax: 0800 90 77 75 Netherlands: Tel: 0800 02 23 042  
Fax: 0800 02 23 043 Germany: Tel: 0800 1 81 00 32 Fax: 0800 1 81 00 31  
United Kingdom: Tel: 0800 89 60 46 Fax: 0800 89 17 14 Switzerland: Tel: 0800 55 50 27  
Fax: 0800 55 50 37 Italy: Tel: 800 79 03 10 Fax: 800 78 03 11 Sweden: Tel: 020 79 07 35  
Fax: 020 79 07 37 Norway: Tel: 800 11 947 Fax: 800 11 948 Denmark: Tel: 80 88 3477  
Fax: 80 88 34 78 Spain: Tel: 900 973 255 Fax: 900 973 256

PC1184-1203