
Tools for SUCCESSFUL

ChIP

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Chromatin was prepared from U2OS cells according to the Abcam X-ChIP protocol. Cells were fixed with formaldehyde for 10 min. The ChIP was performed with 25 µg of chromatin, 2 µg of anti-Histone H3 (tri methyl K4) antibody (ab8580), and 20 µl of Protein A/G sepharose beads. No antibody was added to the beads control (yellow). The immunoprecipitated DNA was quantified by real time PCR (Taaman approach).

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All batches of anti-Histone H3 (acetyl K27) antibody (ab4729) are tested in Peptide Array against peptides to different Histone H3 modifications. Results show strong binding to Histone H3 - acetyl K27 peptide (ab24404), indicating that this antibody specifically recognises the Histone H3 - acetyl K27 modification.

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Product	When to Use	Sample/rxn	Steps covered	Working Time	Cat. No.
ChIP Kit - One Step	Have chromatin ready	5 µg chromatin	From chromatin shearing to recovery of DNA	5 hours	ab117138 (48 / 96 tests)
ChIP Kit Magnetic - One Step	Have chromatin ready	5 µg chromatin	From chromatin shearing to recovery of DNA	5 hours	ab156907 (48 / 96 tests)
High Sensitivity ChIP Kit	Have low input amount	2 x 10 ³ cells 0.5 mg tissue	From chromatin extraction to recovery of DNA	5 hours	ab185913 (24 / 48 tests)
ChIP-Seq High Sensitivity	Want to prepare a NGS library direct from cells	5 x 10 ⁵ cells	From chromatin extraction to recovery of DNA	7 hours	ab185908 (12 / 24 tests)
ChIP Kit - Plants	Have plant tissue samples	20 mg tissue	From chromatin extraction to recovery of DNA	6 hours	ab117137 (24 / 48 tests)

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ChIP performed on HeLa cells using anti-Histone H3 (tri methyl K4) antibody (ab8580) together with ChIP Kit - One Step (ab117138). No antibody was added to the blank (yellow). The immunoprecipitated DNA was quantified by real time PCR (Taqman approach).

High-Sensitivity ChIP Kit

High Abundance Protein Enrichment ab185913

Cell Number	RNA polymerase II antibody	Non-immune IgG
0	1	1
2,000	~10	~1
4,000	~20	~1
20,000	~80	~1
100,000	~300	~1
500,000	~1,000	~1

Sheared chromatin isolated from different numbers of MBD-231 cells was used for ChIP-qPCR analysis of RNA polymerase II enrichment in GAPDH promoters using High-Sensitivity ChIP Kit (ab185913) and a quantitative PCR Fast Kit.

High-Sensitivity ChIP Kit

Low Abundance Protein Enrichment ab185913

Cell Number	ER-alpha antibody	Non-immune IgG
0	1	1
50,000	~25	~1
500,000	~75	~1

Sheared chromatin isolated from different numbers of MCF7 cells was used for ChIP-qPCR analysis of ER-alpha enrichment in TFF1 promoters using High-Sensitivity ChIP Kit (ab185913) and a quantitative PCR Fast Kit.