



abcam



Knockout Cell Lines

Product Catalog



產品諮詢請洽 Abcam 台灣代理 — 伯森生技

歡迎前往 Abcam Knockout Cell Lines 產品專區瀏覽詳細說明與常見問答

abcam.com/reagents/knockout-cell-lines



 伯森生技

網址 blossombio.com 台北 (02) 2788-2121 新竹 (037) 581-106
客服 0800-059668 台中 (04) 2323-3939 南部 (06) 235-6628



Knockout Cell Lines

目錄 Contents

產品特點


- 使用 CRISPR-Cas9 技術剔除單一特定基因
- 採用 Sanger 定序法嚴格驗證把關剔除效果
(諸多商品另附 Western blot 確校實驗結果)
- 超過 450 個產品品項，目前提供 HeLa、HEK293T 兩種細胞株選擇
- 可依據需求免費提供 Wild-type 細胞株
(請於下單時告知伯森生技業務專員)

Homozygous: 1 bp deletion in exon4

```
Mut AGTTCTCTTCATTGACCAAGGAA-TCGGCCTCTATTTGAAGATATGA
WT  AGTTCTCTTCATTGACCAAGGAAATCGGCCTCTATTTGAAGATATGA
```

Human IL18 knockout HeLa cell line (ab265274) is validated by Sanger sequencing.

Mutation: 1 bp deletion in exon4. **Editing tool:** CRISPR/Cas9.
Cell line: Diploid of HeLa cell.



Human IL18 knockout HeLa cell line (ab265274) is validated by Western blot.

Green: Anti-IL-18 antibody [EPR19954-188] (ab243091) observed at 22 kDa.

Red: Loading Control, Anti-GAPDH antibody [6C5] - (ab8245) observed at 37 kDa.

依剔除基因名稱排序

(點擊頁碼即可前往該區產品頁面)

A	3
B	5
C	6
D	9
E	10
F	11
G	12
H	13
I	14
J, K	15
L	16
M	17
N	18
O, P	19
R	21
S	22
T	24
U-Z	26

Wild type 27
cell lines

A

Cat. No.	Product	Pack Size
ab265256	Human ABCC1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264665	Human ABCC3 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264730	Human ABHD6 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265612	Human ABL1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264929	Human ABL2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265251	Human ACSL4 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265024	Human ACVR1B knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264911	Human ADAM17 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264709	Human ADCY3 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab266680	Human ADIPOR2 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab261831	Human ADRB1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265481	Human AGO1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265331	Human AK4 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab266182	Human AKAP1 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab265687	Human AKT3 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261763	Human ALPP knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264728	Human AMH knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265077	Human ANTXR1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264700	Human ANXA3 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab267294	Human APBB1 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab266016	Human APH1A knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264732	Human APLNR knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255362	Human APP knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab265118	Human APTX knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264939	Human ARF1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab267258	Human ARFIP1 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab265797	Human ARHGEF2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261802	Human ARL13B knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264701	Human ARNTL knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab266116	Human ARRB2 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml

A

Cat. No.	Product	Pack Size
ab264699	Human ASIC1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264989	Human ASS1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264908	Human ATF3 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab266457	Human ATF4 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab261800	Human ATF6 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264684	Human ATG14 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261772	Human ATG16L1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261773	Human ATG16L1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264720	Human ATP2B4 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261810	Human AXL knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265392	Human AXL knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml

B

Cat. No.	Product	Pack Size
ab264843	Human BAD knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265277	Human BAK1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255363	Human BAX knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264925	Human BAZ1A knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261832	Human BBC3 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261797	Human BCL10 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264724	Human BCL3 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265410	Human BCL6 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265659	Human BECN1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261767	Human BMP2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264666	Human BMP2K knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264776	Human BMPR1A knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264696	Human BMPR2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265373	Human BRAF knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265002	Human BRD7 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264928	Human BRE knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264705	Human BRSK1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265527	Human BST2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml

C

Cat. No.	Product	Pack Size
ab265143	Human CAB39 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255368	Human CANX knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab265310	Human CASP1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255370	Human CASP3 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab266853	Human CASP8 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261816	Human CAT knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255371	Human CAV1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264711	Human CBLB knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264950	Human CBS knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265439	Human CCDC137 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255372	Human CCL2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265379	Human CCL20 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255348	Human CCND1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265471	Human CCR7 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255373	Human CD14 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255344	Human CD274 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261739	Human CD28 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255374	Human CD36 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265951	Human CD38 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab266324	Human CD47 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab265947	Human CD58 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab256498	Human CD68 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255375	Human CD9 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab267254	Human CDC16 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab264731	Human CDC25A knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab266522	Human CDC42 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab255376	Human CDH1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255377	Human CDH2 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab264672	Human CDK16 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264985	Human CDK17 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml

C

Cat. No.	Product	Pack Size
ab264668	Human CDK19 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab255378	Human CDK4 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab266059	Human CDK6 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab265087	Human CDK8 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab255349	Human CDKN1A knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab264877	Human CDKN1B knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab264906	Human CDKN1C knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab255351	Human CDKN2A knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab261761	Human CDKN2B knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab265031	Human CDKN2C knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab261771	Human CEBPB knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab266018	Human CHAT knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab264815	Human CHEK2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab265053	Human CHUK knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab264678	Human CIC knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab265464	Human CLCN2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab267270	Human CLK2 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1 ml
ab264693	Human CLK3 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab265655	Human CLPX knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab265620	Human CNN1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab264949	Human CNP knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab265142	Human CNR1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab261757	Human COIL knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab256349	Human COL1A1 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1 ml
ab255379	Human COL4A1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab264707	Human COL9A3 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab266380	Human COPS7A knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1 ml
ab261809	Human CPS1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab266319	Human CPT1A knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1 ml
ab265931	Human CPT2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml

C

Cat. No.	Product	Pack Size
ab264752	Human CREB5 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255381	Human CRP knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265786	Human CSAD knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265273	Human CSF1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab267226	Human CSF1R knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264667	Human CSNK1D knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264674	Human CSNK1E knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265235	Human CSRP1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255382	Human CTGF knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265893	Human CTLA4 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255352	Human CTNNB1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab266521	Human CTSL knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab261774	Human CXCL1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255383	Human CXCR4 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml

D

Cat. No.	Product	Pack Size
ab266263	Human DAG1 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab266008	Human DCAF13 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265115	Human DDB2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265760	Human DDIT3 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265915	Human DDR1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264864	Human DDR2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261765	Human DKK1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab267290	Human DLG1 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab255384	Human DMD knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261793	Human DNMT3A knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264771	Human DPP4 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265949	Human DTX2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265447	Human DUSP6 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265198	Human DYRK3 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264764	Human DYSF knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml

E

Cat. No.	Product	Pack Size
ab265362	Human E2F3 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264738	Human ECE2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265355	Human EDN1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261822	Human EEA1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265232	Human EFNA1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261830	Human EGF knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255385	Human EGFR knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265149	Human EHMT2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab266208	Human ELOF1 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab265178	Human ENG knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255386	Human EPCAM knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264760	Human EPHA4 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab266732	Human EPHB4 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab255387	Human ERBB2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264881	Human ERCC6 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255388	Human ERG knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255389	Human ERN1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab266543	Human ERP29 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab264714	Human ETV5 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264912	Human EZR knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml

F

Cat. No.	Product	Pack Size
ab261817	Human FADD knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264723	Human FAM3C knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261743	Human FANCD2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255353	Human FGA knockout A549 cell line	1 x 10 ⁶ cells/vial, 1ml
ab265974	Human FGF21 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265268	Human FLCN knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255390	Human FLT1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab267227	Human FLT3 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255391	Human FN1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264921	Human FOLR1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261823	Human FOXA1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265456	Human FOXO4 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265533	Human FTL knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265534	Human FTL knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab266587	Human FUS knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab266133	Human FYN knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab265463	Human FZD1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264727	Human FZD2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml

G

Cat. No.	Product	Pack Size
ab266259	Human GAA knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab264992	Human GAL knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265027	Human GALT knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265038	Human GBA knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255393	Human GCG knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264695	Human GDA knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab256501	Human GFAP knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265560	Human GLOD4 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264968	Human GNAI1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264689	Human GNB1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab266505	Human GOLM1 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab264715	Human GPSM2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264769	Human GRIN1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264811	Human GRK6 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265201	Human GSN knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264703	Human GSTT2B knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml

H

Cat. No.	Product	Pack Size
ab256499	Human H3F3A knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264716	Human HAO1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265120	Human HAVCR2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264690	Human HDAC5 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264804	Human HDAC6 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265221	Human HGF knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255394	Human HIF1A knockout HCT116 cell line	1 x 10 ⁶ cells/vial, 1ml
ab265195	Human HIPK2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261794	Human HIST1H1C KO HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab267279	Human HK1 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab255395	Human HMGB1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab266358	Human HMGB2 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab255396	Human HMOX1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab266657	Human HNRNPA1L2 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab266404	Human HNRNPA2B1 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab265868	Human HPDL knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab266313	Human HSP90B1 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab261738	Human HSPB1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264710	Human HTR2A knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265976	Human HTT knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml

Cat. No.	Product	Pack Size
ab261755	Human IDE knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264916	Human IDH1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264747	Human IDS knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265698	Human IFITM3 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265274	Human IL18 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265871	Human IL1R1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265767	Human IL6ST knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255397	Human INS knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255345	Human IRF3 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265672	Human ITGA5 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265838	Human ITGA6 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255354	Human ITGB1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265353	Human ITGB3 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265847	Human ITPR1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml

J, K

Cat. No.	Product	Pack Size
ab265424	Human JAG1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265864	Human JAK2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab266272	Human JUP knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab264735	Human KCNJ11 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265790	Human KDM1A knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265110	Human KDM6A knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab266294	Human KDM6B knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab264742	Human KIF7 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264768	Human KIT knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261791	Human KLF5 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab267229	Human KRAS knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265010	Human KRT10 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255398	Human KRT14 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255399	Human KRT18 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265706	Human KRT20 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255400	Human KRT8 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml

L

Cat. No.	Product	Pack Size
ab255401	Human L1CAM knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261798	Human LAG3 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255402	Human LAMP2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264719	Human LAMP3 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255403	Human LDHB knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab265327	Human LDLR knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261780	Human LGR5 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264737	Human LIF knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261787	Human LMNA knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255404	Human LMNB1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261801	Human LOX knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261804	Human LOXL2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264807	Human LOXL2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265482	Human LPL knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab266618	Human LRP5 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab264726	Human LRRK1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255405	Human LRRK2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265153	Human LTB4R2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml

Cat. No.	Product	Pack Size
ab264743	Human MAP3K12 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab264669	Human MAP4K4 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab255406	Human MAPK14 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1 ml
ab265042	Human MAPKAPK3 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab266653	Human MAPRE1 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1 ml
ab264722	Human MC2R knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab261790	Human MCAM knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab265961	Human MET knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab265821	Human MGEA5 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab255407	Human MKI67 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab264758	Human MKNK1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab267223	Human MLH1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab255408	Human MLKL knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab264794	Human MME knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab265872	Human MMP3 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab255409	Human MMP9 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab264704	Human MOK knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab264679	Human MPI knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab264878	Human MPLKIP knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab265715	Human MPP5 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab255410	Human MSH6 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab261828	Human MSLN knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab265272	Human MTAP knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab255411	Human MTOR knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1 ml
ab255412	Human MUC1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab256500	Human MYC knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1 ml
ab266036	Human MYL9 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1 ml
ab267297	Human MYSM1 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1 ml

N

Cat. No.	Product	Pack Size
ab264756	Human NBL1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261834	Human NBN knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265275	Human NCAM1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255413	Human NEFH knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264755	Human NEK1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264856	Human NEK4 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264725	Human NF1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265199	Human NFATC1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255414	Human NFE2L2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264770	Human NGF knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265001	Human NLK knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264955	Human NME2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab266244	Human NONO knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab255416	Human NOS2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265099	Human NPHS1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261814	Human NPPA knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261795	Human NPY knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261766	Human NR3C1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265540	Human NR3C2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265334	Human NR4A1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264734	Human NRG4 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261808	Human NTRK2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264759	Human NUPR1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml

O, P

Cat. No.	Product	Pack Size
ab264749	Human ORAI3 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264713	Human OSMR knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265271	Human PARK2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264973	Human PCNT knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264775	Human PCSK9 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255417	Human PDCD1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265191	Human PDCD1LG2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265628	Human PDCD1LG2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261833	Human PDCD4 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265342	Human PDE10A knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264729	Human PDE3A knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264718	Human PDE3B knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264740	Human PDGFC knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265304	Human PDGFD knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264791	Human PDK1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261805	Human PDK4 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab256502	Human PECAM1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264741	Human PFKFB4 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab266766	Human PGM1 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab261785	Human PGR knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264902	Human PHF6 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264683	Human PI4K2B knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264691	Human PKN2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265021	Human PKN3 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265255	Human PLAA knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261811	Human PML knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261776	Human PMS2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264984	Human PODXL knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261748	Human POLB knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264754	Human POLL knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml

O, P

Cat. No.	Product	Pack Size
ab261806	Human POSTN knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab266542	Human POU4F2 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab255418	Human POU5F1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265631	Human PPARD knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261812	Human PPARG knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261821	Human PPARG knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab266162	Human PPM1G knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab264673	Human PRKACB knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265721	Human PRKCD knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264670	Human PRKCI knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264839	Human PRMT3 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264676	Human PSCA knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255419	Human PTEN knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255420	Human PTGS2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255421	Human PTK2 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab255422	Human PTPRC knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab266231	Human PTTG1 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab265285	Human PUM1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264767	Human PXN knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml

R

Cat. No.	Product	Pack Size
ab255423	Human RAB7A knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264993	Human RAB8A knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264994	Human RAB8A knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264978	Human RAF1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265092	Human RALA knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265404	Human RALBP1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264677	Human RBBP7 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255424	Human RBFOX3 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264706	Human RBL1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab266663	Human RBM3 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab264697	Human RBPMS knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265478	Human RCHY1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab266975	Human REC8 knockout A549 cell line	1 x 10 ⁶ cells/vial, 1ml
ab255425	Human RELA knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab266154	Human RHOC knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab264688	Human RIPK2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264761	Human RRAGB knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264753	Human RRAGC knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261783	Human RUNX2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml

S

Cat. No.	Product	Pack Size
ab261758	Human S100A4 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255426	Human S100B knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab267038	Human SAMD9 knockout A549 cell line	1 x 10 ⁶ cells/vial, 1ml
ab265567	Human SCAMP1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264733	Human SCN1B knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264750	Human SERPINB5 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261779	Human SFN knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264746	Human SIAH1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265933	Human SIKE1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265434	Human SIPA1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255428	Human SIRT1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265133	Human SLC12A6 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265668	Human SLC25A13 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261747	Human SLC2A3 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265262	Human SLC2A4 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265674	Human SLC2A8 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264664	Human SLC30A1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264762	Human SLC35B2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264686	Human SLC39A8 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab266234	Human SLC7A1 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab255430	Human SMAD2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255431	Human SMAD3 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264745	Human SMAD7 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255432	Human SMARCA4 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab264859	Human SMARCC1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265963	Human SNAI1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255433	Human SNCA knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab255434	Human SOD2 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab264757	Human SP100 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265741	Human SPG7 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml

S

Cat. No.	Product	Pack Size
ab265742	Human SPG7 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255435	Human SPP1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264681	Human SPR knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265787	Human SQLE knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255343	Human SQSTM1 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab265016	Human SRC knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265911	Human SRGN knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265443	Human SRPK1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264865	Human SRPK2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265151	Human SSTR1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265236	Human SSTR2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255346	Human STAT1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261819	Human STAT2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255436	Human STAT3 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265640	Human STAT5A knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265442	Human STK4 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab266005	Human SYNJ1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265359	Human SYNJ2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255356	Human Syp knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab267272	Human Syp knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml

T

Cat. No.	Product	Pack Size
ab266845	Human TAB1 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab264717	Human TBL1X knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264828	Human TBXA2R knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265127	Human TCF7 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab267064	Human TCHH knockout A549 cell line	1 x 10 ⁶ cells/vial, 1ml
ab265437	Human TET2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264946	Human TFE3 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255438	Human TFRC knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255439	Human TGFB1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265245	Human TGM2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264880	Human THRA knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255440	Human THY1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265059	Human TLE1 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab264901	Human TLE1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264898	Human TLK1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265743	Human TM4SF1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265227	Human TMEM173 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264744	Human TMPRSS6 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265852	Human TMUB1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265511	Human TNC knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265983	Human TNFAIP3 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264687	Human TNFRSF10A knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265856	Human TNFRSF11B knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265347	Human TNFRSF18 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265460	Human TNFRSF8 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab261825	Human TNFSF10 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264748	Human TOB1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255441	Human TOMM20 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab261759	Human TP53 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265406	Human TP63 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml

T



Cat. No.	Product	Pack Size
ab264963	Human TRIM24 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255442	Human TRIM28 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265079	Human TRRAP knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264698	Human TSC22D4 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265901	Human TSHR knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264942	Human TSPO knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255358	Human TUBB3 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265138	Human TYMP knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml

U-Z

Cat. No.	Product	Pack Size
ab264842	Human UBE2D1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264777	Human UBE2D2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab266346	Human UBE2E1 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab255443	Human UCHL1 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab264692	Human UHMK1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265725	Human UNG knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265015	Human USP2 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264682	Human USP48 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab266293	Human VAMP8 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab265580	Human VCL knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264680	Human VCIPI1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255444	Human VDAC1 knockout HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab255445	Human VEGFA knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255446	Human VIM knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264721	Human VRK3 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265483	Human VTI1B knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264663	Human XPA knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264739	Human XRCC4 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265202	Human YES1 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab265522	Human ZAP70 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab266951	Human ZC3HAV1 knockout A549 cell line	1 x 10 ⁶ cells/vial, 1ml
ab264736	Human ZMPSTE24 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264671	Human ZMYM3 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264874	Human ZMYND8 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab264694	Human ZNF622 knockout HeLa cell line	1 x 10 ⁶ cells/vial, 1ml

Wild Type Cell Lines

Cat. No.	Product	Pack Size
ab255450	Human Wild Type A549 cell line	1 x 10 ⁶ cells/vial, 1ml
ab255451	Human Wild Type HCT116 cell line	1 x 10 ⁶ cells/vial, 1ml
ab255449	Human Wild Type HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab255593	Human Wild Type HEK293T cell line	1 x 10 ⁶ cells/vial, 1ml
ab255448	Human Wild Type HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab255928	Human Wild Type HeLa cell line	1 x 10 ⁶ cells/vial, 1ml
ab260075	Human Wild Type HeLa cell line	1 x 10 ⁶ cells/vial, 1ml



產品諮詢請洽 Abcam 台灣代理 — 伯森生技

歡迎前往 Abcam Knockout Cell Lines 產品專區瀏覽詳細說明與常見問答

abcam.com/reagents/knockout-cell-lines



伯森生技

網址 blossombio.com 台北 (02) 2788-2121 新竹 (037) 581-106
客服 0800-059668 台中 (04) 2323-3939 南部 (06) 235-6628

