


Western Blot Stripping Buffer

Cat. No.: ab270550 | Pack Size: 500 ml, 5 L

- 不需加熱步驟，室溫 20 分鐘即可迅速完成抗體脫附！
- 效果穩定，不會造成轉漬膜上的抗原損傷。
- 可進行多次抗體探測，有效節省珍貴樣本。




Saves sample by re-probing. Reprobing of the same membrane with different antibodies after treatment with Western Blot Stripping Buffer (ab270550).

Bradford Reagent

Cat. No.: ab119216 | Pack Size: 500 ml

- 專利改良配方，定量結果更準確，不易受到樣本中的 Detergent 影響！
※高含量定量流程，可耐受 1% Detergent；微量定量流程，可耐受 0.1% Detergent。
- 偵測範圍：0.1~1.5 mg/ml (高含量定量) 或 1~25 µg/ml (微量定量)。



Abcam Bradford Reagent is not affected by the detergents compared to leading competitor "Coomassie Protein Assay" which are subject to the high background noise at 595 nm in the presence of detergents.

Western Blot Stripping Buffer 與 Bradford Reagent 皆為 Expecon (part of Abcam) 所開發的高效蛋白質實驗分析試劑。完整產品資訊與最新活動訊息，請洽 Abcam 台灣代理 — 伯森生技。



相關文獻：

Bradford Reagent (ab119216)
previously called Bradford Ultra (BFU05L, Expecon)

- PLRP2 selectively localizes synaptic membrane proteins via acyl-chain remodeling of phospholipids. *J Lipid Res.* 2020 Sep 22; jlr.RA120001087. PMID: 32963038
- Bacterial polyphosphates interfere with the innate host defense to infection. *Nat Commun.* 2020 Aug 12; 11(1):4035. PMID: 32788578
- The RNA-Binding Protein A1CF Regulates Hepatic Fructose and Glycerol Metabolism via Alternative RNA Splicing. *Cell Rep.* 2019 Oct 8; 29(2):283-300.e8. PMID: 31597092
- Aquaporins influence seed dormancy and germination in response to stress. *Plant Cell Environ.* 2019 Aug; 42(8):2325-2339. PMID: 30986891

相關連結：



Western Blot Stripping Buffer (ab270550)



Bradford Reagent (ab119216)

