

abcam

# 螢光一級抗體

## 8折特惠

活動時間  
2022/4/25~6/29

### 豐富齊全的螢光標定一級抗體選擇

- 超過 5,000 種螢光標定一級抗體選擇
- 適用流式細胞分析 (Flow cytometry)、螢光染色 (IF/ICC/IHC)、螢光西方墨點偵測 (WB) 等實驗
- 滿足多目標染色 (multiplex)、共定位 (co-localization) 實驗需求



Abcam Conjugated  
Primary Antibodies

螢光標定分子 Ex/Em		DyLight® 550 562/576	PE/DyLight® 594 593/618	Spectral Red™ 565/666	Alexa Fluor® 680 679/702
Alexa Fluor® 405 402/421	Alexa Fluor® 488 499/519	PE 496, 566/576	SureLight® APC 652/657	DyLight® 650 654/673	redFluor™ 710 640/710
violetFluor™ 450 405/450	Cy3® 548/561	Alexa Fluor® 568 578/603	APC 650/660	PerCP/Cy5.5® 482/676	Alexa Fluor® 700 702/723
CF405M 408/452	Alexa Fluor® 555 555/565	Texas Red 595/613	Cy5® 647/665	PerCP 477/678	APC/Cy7 650/774
DyLight® 488 495/519	Rhodamine 550/570	PE/Texas Red® 566/616	Alexa Fluor® 647 650/665	PE-Cy5.5® 565/693	PE-Cy7 566/778
FITC 495/519	TRITC 557/576	Alexa Fluor® 594 590/617	PE-Cy5® 565/666	APC-Cy5.5® 650/694	Alexa Fluor® 790 782/805

產品資訊 · 活動諮詢

請洽 Abcam 台灣代理伯森生技



Multicolor  
flow cytometry  
實驗設計技巧與  
注意事項



Multicolor  
flow cytometry  
抗體組合快選工具

伯森生技

伯森生技客服專線

北部 0800-059668 南部 0800-259988

伯森網站



FB



LINE

